Fig. 1


Fig. 2


Fig. 3



Fig. 4

Product of Mexico. 8362765, D12556-3 B(7-10)

## EN SUGGESTIONS FOR USE

### CLINICIANS MUST BE TRAINED IN THE PRACTICE OF VENIPUNCTURE AND BE AWARE OF INHERENT DANGERS.

Aseptic technique, proper skin preparation and continued protection of the site are essential. Observe Universal Precautions on **ALL** patients.

**SLIDE CLAMP SHOULD NOT BE ENGAGED UNTIL WIRE STYLET IS WITHDRAWN.**

- A. Remove needle cover and inspect unit. Secure flow control plug or remove to preprime/preconnect system.
- B. Release catheter seal by securing wing with fingers and pushing upward on the segment directly below wings (fig. 1) Check bevel orientation and lie distance.
- C. Grasp textured side of wings, pinching firmly, and make venipuncture (fig.2).
- D. Observe blood return and advance catheter into vein.

### NEVER REINSERT THE NEEDLE INTO THE CATHETER. DO NOT BEND THE NEEDLE OR STYLET WHILE USING THE PRODUCT.

- E. Apply pressure to middle of each wing. Grasp textured end of shield activator and withdraw needle stylet in a straight continuous motion (fig. 3.) **Important: Do not use force to withdraw needle stylet. If necessary remove entire unit and discard.**
  - F. Discard shielded needle stylet immediately in a puncture resistant, leak proof sharps container.
  - G. Connect tubing or adapter. **Flow control plug should always be replaced and never used as a cap.**
  - NOTE:** If shield has not locked over needle tip, grasp textured end of shield activator with needle tip pointed down. Gently push downward until needle is covered (fig. 4). **Keep needle point away from body and fingers at all times.**
  - H. Secure catheter using H-method and dress site per protocol.
- REPORT NEEDLESTICKS IMMEDIATELY AND FOLLOW ESTABLISHED PROTOCOL.**
- Percutaneous puncture with a contaminated needle may lead to serious illness such as hepatitis, HIV (AIDS) or other infectious diseases.
  - Resheathing needles is hazardous.

**NON PYROGENIC-STERILE:** Unless package has been opened or damaged. **DO NOT** use scissors at or near insertion site.

This product does not contain natural rubber latex. Re-use may lead to infection or other illness/injury.

**CAUTION:** Use of this product with power injectors may cause catheter leakage or product damage.

This product contains a phthalate which may cause a risk to {children (male), pregnant or nursing women} in conjunction with prolonged duration of use.

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|                                                                                       |                                                     |
|---------------------------------------------------------------------------------------|-----------------------------------------------------|
| <b>STERILE</b> <b>EO</b>                                                              | Sterilized using ethylene oxide                     |
| <b>LOT</b>                                                                            | Batch code                                          |
|    | Caution, consult accompanying documents             |
|    | Do not reuse                                        |
|    | Use by                                              |
|    | Green dot                                           |
|    | Recyclable packaging                                |
|   | Recyclable packaging                                |
| <b>REF</b>                                                                            | Catalogue number                                    |
|  | Do not use if package is damaged                    |
|  | Medical device containing phthalates                |
|  | Manufacturer                                        |
| <b>EC REP</b>                                                                         | Authorised representative in the european community |

BD and BD Logo are trademarks of Becton, Dickinson and Company. © 2010 BD.

|                                                                                       |                                                                                                 |
|---------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------|
|  | Becton Dickinson Infusion Therapy Systems Inc., 9450 South State Street, Sandy, Utah 84070, USA |
| <b>EC REP</b>                                                                         | BD Distribution Centre, Laagstraat 57, B-9140 Temse, Belgium                                    |

Becton Dickinson Korea Ltd., 5th & 6th Floor, Im Sung Building, 788-2, Yoksam-Dong, Kangnam-Ku, Seoul 135-080, Republic of Korea
Becton Dickinson India Private Limited, The Surya, Best Western Hotel, E Floor, New Friends Colony, New Delhi, 110 065, India
Becton Dickinson de Mexico, S.A. de C.V., Monte Pelvoux ill 8° Piso, Lomas de Chapultepec, 11000 Mexico, D.F.
Becton Dickinson Philippines, Inc., 29F, The Discovery Center, #25 ADB Avenue, Ortigas Center, Pasig City
Becton Dickinson & Company, Singapore Branch, 30 Tuas Avenue 2, Singapore, 639461
Becton Dickinson Canada Inc., 2771 Bristol Circle, Oakville, Ontario, Canada L6H 6R5
Becton Dickinson Pty Ltd, 4 Research Park Drive, North Ryde NSW 2113 Australia
Becton Dickinson Ltd 8 Pacific Rise Mt Wellington, Auckland New Zealand
Becton Dickinson Ind. Cir. Ltda - Juiz de Fora - MG – Brazil
福島県福島市土船字五反田 1 番地

**BD Medical**

9450 South State Street  
Sandy, Utah 84070 USA.

For information not listed here, see appropriate Buy Specification.

**Title:** BD Saf-T-Intima Multiuse IFU

**Artwork Drawing Number:** D12556

**Single Unit Dimensions:** 11" x 17"

**Print Colors:** PMS 2755 Blue

**Notes:** Finished folded size is 5.5" x 5.67" with the BD Logo panel facing up.

**Material:** 50# White Offset

**Revision History**

| Rev | Description of Changes                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         | ECO No.   |
|-----|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------|
| 2   | <p>Update labeling to meet the requirements found in the new Medical Device Directive (MDD) per PDN 08-3300-338.</p> <p>1.) Add the Recycle symbol "circle arrow".</p> <p>2.) Add green dot "recycle symbol".</p> <p>3.) Add recycle symbol "triangle arrow".</p> <p>4.) Remove latex symbol and moved statement This product does not contain natural rubber latex to main instruction section.</p> <p>5.) Add the Manufacture symbol in front of the Sandy Address.</p> <p>6.) Add the EC rep symbol in front of the Temse Address.</p> <p>7.) Update addresses from S-Drive.</p> <p>8.) Add Rx Only.</p> <p>9.) Add the DEHP symbol.</p> <p>10.) Patent Numbers (need to check with Legal).</p> <p>11.) Add the Risk re-use statement "Re-use may lead to infection or other illness/injury".</p> <p>12.) Add the package damage symbol.</p> <p>13.) Add sterile Eto symbol.</p> <p>14.) Add CE-Mark and Single Use symbols.</p> <p>15.) Add the explanation of all symbols.</p> <p>Update to current Dgp Format and add history page. ECO # will be 96542</p> <p>Material number D12556C will be replaced by SAP # 8359303. Version will go to 2 B(4-09)</p> <p>Removed Japan language. Federal (USA) Law restricts..... needs to be removed from all languages. This will be replaced with the Rx only symbol mentioned in the other list of changes.</p> | ECO96542  |
| 3   | <p>Replace 8359303 with new material number 8362765; update file name to D12556-3 and date to 5-10. Insert correct illustrations for figures 1 and 3. Figure 1 must include the clamp as part of the visual aid and Figure 3 must reflect an STI device instead of an Intima device in the visual aid. Add statement, "This medical device contains a phthalate which may cause a risk to {children (male), pregnant or nursing women} in conjunction with prolonged duration of use.</p> <p>Update trademark by-line.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     | ECO122805 |
|     |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |           |