


# ● Tubigrip™ Shaped Support Bandage (TSSB)

## ANATOMICALLY SHAPED BANDAGE

- TSSB is a shaped support bandage used to aid venous and lymphatic return in the management of venous disorders of the legs and arms

- Washable and reusable
- Easy to apply and re-apply
- Choice of 2 colours for patient acceptability



Full range of sizes to suit most applications

Anatomically shaped to provide support along the complete length of the limb

Provides effective light compression<sup>1</sup>

## HOW TSSB WORKS

Tubigrip Shaped Support Bandage provides lasting, effective support, while allowing ambulant patients complete freedom of movement. The bandage only stretches radially and this property enables it to provide support without rucking or slipping down. TSSB can be washed for repeated use.

Only a single measurement around the calf or forearm is required for precise size selection. The range of five full-leg and three below-knee bandages ensure a good fit for the majority of patients. The choice of two shades helps suit an individual patient's preference – natural and flesh.

## BENEFITS OF TSSB

- Anatomically shaped to provide support along the complete length of the limb
- Elasticated cotton fabric construction is comfortable as well as effective
- Simple size selection is quick and convenient
- Washable and reusable
- Reapplication after washing is quick and easy for patients.

## MATERIAL PROPERTIES

- Cotton
- Elastodiene (contains natural rubber latex)
- Polyamide

## AREAS OF USE

Tubigrip Shaped Support Bandage (TSSB) is a support bandage used to aid venous and lymphatic return in the management of venous disorders of the legs and arms.

TSSB may also be used as an undercast stockinette in cast bracing, for support following cast removal, and to provide pressure therapy in the treatment of post-burn scarring.

### Washing instructions



### TSSB ORDERING INFORMATION

Art. No	Size	Limb cm	Pieces Per Case	PIP Code	NHS Code
<b>Natural shade Below-knee</b>					
1472-03	Small B/C	32-36cm	60	054-5517	EGB011
1473-03	Medium C/D	35-39cm	60	054-5509	EGB013
1474-03	Large D/E	38-42cm	60	054-5491	EGB015
<b>Full-leg</b>					
1475-03	Small B/C	32-36cm	110	054-5483	EGB031
1476-03	Medium C/E	35-39cm	110	054-5475	EGB500
1477-03	Large D/F	38-42cm	110	054-5467	EGB035
1478-03	X-Large E/G	38-42cm (For thighs over 60cm)	110	054-5459	EGB037
1471-03	X-Large X/T	41-45cm (For thighs Over 60cm)	110	074-2841	-
<b>Flesh Below-Knee</b>					
1479-03	Small B/C	32-36cm	60	076-6337	-
1480-03	Medium C/D	35-39cm	60	076-6329	-
1481-03	Large D/E	38-42cm	60	076-6311	-
<b>Full-leg</b>					
1482-03	Small B/D	32-36cm	110	076-6303	-
1483-03	Medium C/E	35-39cm	50	076-6295	-
1484-03	Large D/F	38-42cm	50	076-6287	-
1485-03	X-Large E/G	38-42cm (For thighs over 60cm)	110	076-6279	-
1486-03	X-Large X/T	41-45cm (For thighs over 60cm)	110	309-4117	-

Individually packed incartons of 10.

The below knee style is also suitable for use on the arms.

### References:

1. Cornwall JV, Dore CJ, Lewis JD. Graduated compression and its relation to venous refilling time. British Medical Journal 1987; 295: 1087-1090



**MÖLNLYCKE  
HEALTH CARE**

Mölnlycke Health Care Ltd, Arenson Centre, Arenson Way, Dunstable,  
Bedfordshire LU5 5UL, UK. Tel. +44 (0) 870 6060766, [www.molnlycke.co.uk](http://www.molnlycke.co.uk)

The Mölnlycke Health Care name and logo, Tubigrip™ are registered trademarks of Mölnlycke Health Care AB.

©2012 Mölnlycke Health Care. All rights reserved.

UK162550611

