


QUICK USER GUIDE

Contec Infrared Thermometer



Code: 27125


Step 1



Press the “on” button once quickly and release pressure off the trigger



Step 2



Ensure the screen shows that the unit is in Body mode

Body mode

Step 3

Position the unit vertically aligning the beam to the centre of the forehead (between the eyebrows). It must be less than 3cm's from the patient - do not touch the patient.

Step 4

Press the “on” button once quickly then release the pressure on the trigger.



Step 5



View the measurement on the screen

The screen backlight will change in colour to indicate temperature rating status as below:



Mild temperature shown in yellow

Within normal range shown in blue

High temperature shown in red

Review memory for past readings



From the off state press and hold trigger button to view first 4 readings

Click once to scroll through to other saved screens memory capacity is 30 groups of data

To enter settings mode



Repeat as for memory setting once in memory mode hold down trigger and setting mode will be displayed.

User consideration:

- Before measuring ensure the measurement area (forehead) is not covered by hair, sweat, moisture or cosmetics as this can interfere with measurement outcomes
- Forehead temperature measurements can be influenced by environmental temperature. If the forehead is sweaty or clammy measure by aiming the unit behind the earlobe
- Ensure 15 seconds between each measurement
- The ambient temperature around the unit should be stable, do not measure in places with large air flow, such as fan air conditioning outlets
- Ensure the thermometer is in the correct mode which should always be Body mode and is visible on the measurement screen. Obj mode is a calibration mode and will not work for taking body temperatures
- Keep the thermometer a 1 metre distance away from electromagnetic items eg: mobile phones when performing a measurement