



- **OVERCOME THE PROBLEMS** of high cost, cross infection risks and health and safety hazards commonly associated with reprocessing reusable LLETZ electrodes.
- **Rocket® LLETZ LOOPS** are designed to deliver higher cutting performance with less thermal damage to the sample.
- **2.4mm (3/32") x 10cm shaft** to fit all standard UK, US and European diathermy systems
- **DELTA™** style loops available for increased wire rigidity and increased control of penetration.

Rocket® LLETZ Loops were the world's first single use loop electrodes designed specially for LLETZ.

Rocket® LLETZ Loops have been in continuous production for more than 27 years.

They continue our tradition of providing exceptional value, reliability and quality

Rocket® LLETZ Range	Size		Code
Rocket® LLETZ Diathermy Loop, for single use, packed sterile in cartons 25 units	Small	14mm	R57800
	Medium	18mm	R57801
	Large	22mm	R57802
	Square	10mm	R57803
Rocket® DELTA™ LLETZ Diathermy Loop, for single use, packed sterile in cartons 25 units	Medium	18mm	R57851
	Large	22mm	R57852
Rocket® LLETZ Ball Electrodes, for single use, packed sterile in cartons 25 units	Standard	5mm	R57804
	Small	3mm	R57807
Rocket® LLETZ Diathermy Extension, for single use, individually packed sterile		50mm	R57818
Rocket® NETZ™ Needles, for single use, packed sterile in cartons 25 units		15mm	R57860