CLEAN SLEEP PILLOW

This Pillow is covered in 'Clean Sleep' fabric, which enables the product to be

- Waterproof
- Anti-bacterial

Flame Retardant

Plus the added feature of the "Mardel Valve", which enables the pillow to breathe and ensure patient comfort.



Microfresh

Microfresh is a revolutionary product which has been applied to the fibre fill of Pillows. It works to kill a broad spectrum of bacteria including Staphylococcus aureus, fungi, moulds and house dust mites.

Microfresh helps prevent the spread of disease and infection caused by contact with items that have been contaminated by sneezing or coughing. Microfresh also maintains odour-free surfaces by killing the moulds which produce the smells.



Physical Fabric Properties Specification

ТҮРЕ	TEST METHOD	RESULT
Tear Strength, (N/mm)	AS 2324 – 1979	17.5 N/mm minimum
Dimensional Stability @ 80C, (%)	AS 2324 – 1979	4.0% maximum
Hardness, BSS	BS 2782	45 +/- 5
Elongation, (%)	AS 2324	200% minimum
Tensile Strength, MPa	AS 2324	11.0 MPa minimum
Flammability	BS 1763	10mm maximum
Bacterial Inhibition	AATCC - 90 - 1970	O, O HALO (mm) Satisfactory

Cleaning Recommendations

- Removal of typical soiling with warm soapy water or very mild vinyl cleaner is recommended.
- Launder and tumble drying are not recommended.
- Phenolic disinfectants, abrasive non-aqueous product or dry-cleaning solvents (eg: trichloroethane types) are not to be used.
- Disinfectants of a non-phenolic composition are recommended in a diluted solution.

