

Atrauman[®] Silicone

The new silicone contact layer
from HARTMANN



Soft, can be cut to size,
thin, can be left in place
up to 7 days.

The very thin weave of Atrauman Silicone prevents the newly formed granulation tissue growing through the dressing, while allowing the easy passage of exudate through to the secondary dressing.

A wide range of indications

Skin Tears



Burns



Abrasions



**Grafts &
Donor sites**



NPWT



Challenging areas



Provides secure, atraumatic protection

The new Atrauman Silicone is specifically designed for sensitive skin, to provide secure and atraumatic wound care.



Benefits

- Maximum patient comfort thanks to atraumatic and pain-free dressing changes.
- Easy to use: can be cut, very thin and conformable for optimal contact with the wound bed.
- Cost effective: can be left in place for up to 7 days, often only the secondary dressing needs to be changed.

“Sometimes what appears to be a minor injury can prevent you doing the things you love. Having the right dressing for me has made such a difference.”

-Dave Miller, 72

THREE SIZES



APPLICATION:



Suggested fixation:
Peha-haft, cohesive bandage
Omnifix, non woven tape
Hydrofilm, Waterproof film

Atrauman[®]

A wide range of sizes to suit all your needs

For secure, atraumatic protection:

Atrauman[®] Silicone



Product code	Size	Pack content
17090	7.5 × 10 cm	5 pcs
17091	10 × 20 cm	5 pcs
17092	20 × 30 cm	5 pcs



ORDER BY FREEPHONE:
Auckland: 0800 658 814
Christchurch: 0800 804 546
ORDER BY FREEFAX:
Auckland: 0800 830 660
Christchurch: 0800 804 547

AUCKLAND HEAD OFFICE
494 Rosebank Road Avondale
AUCKLAND NUMBERS:
Telephone: 09 829 0960
Facsimile: 09 829 0962
CHRISTCHURCH NUMBERS:

Telephone: 03 348 5260
Facsimile: 03 348 9830