

WelchAllyn

Valley Internal Medicine 09:37 12/31/2012

NIBP AVERAGED
120/80
SYS/DIA mmHg (MAP 93)
AVERAGING : Program 1

PULSE RATE
73
♥/MIN
SOURCE : SpO2

BMI **20.6** **WEIGHT** **61.3** kg
HEIGHT **172.7** cm

SpO2 **97%** (PI 19.3)

TEMPERATURE **38.6** °C (101.5° F)


PATIENT **David A. Barker** **Adult** **PAIN** **2** **Clear** **Save**

Home **Review** **Settings**

WelchAllyn
25 Thermometer
Probe Covers

Manufactured by Welch Allyn, Inc.
Branford, Connecticut, U.S.A.





Welch Allyn Connex[®] Vital Signs Monitor

Specifications and Ordering Information

FEATURES:

Industry-Leading Parameters

- > SureTemp[®] Plus oral, pediatric axillary, adult axillary and rectal thermometry
- > Braun ThermoScan[®] PRO 4000 ear thermometry
- > Blood pressure averaging
- > Masimo[®] or Nellcor[®] pulse oximetry (SpO₂) with percentage or waveform display
- > Fast 15-second SureBP[®]—takes blood pressure measurements upon inflation for the fastest speed on the market (about 15 seconds) contributing to greater patient comfort
- > Height, weight and BMI measurements
- > Spot check of NIBP, pulse rate, temperature and SpO₂ in approximately 30 seconds
- > Available in mobile, wall-mounted or desktop workflow options

USB Accessory Ports

Depending on the model you choose, one or four USB ports are included with your monitor to allow for additional accessory plug-ins such as the accessory cable management stand LED light, the Braun ThermoScan PRO 4000 Ear Thermometer, or select Health o meter[®] Digital Weight Scales.

Partners in Care ServicesSM

Protection throughout the life of your investment, with flexible programs that cover all of your service needs. While product warranties provide basic assurance of Welch Allyn hardware quality, they may not include the full range of services and support you need to achieve your service level and cost objectives. Comprehensive Partnership Programs are available in one- and two-year service agreements.

For more information on Partners in Care Services please visit www.welchallyn.com/services or ask your Welch Allyn representative for more details.



Configurations available:

Connex VSM	6300	6400
Parameters		
SureBP NIBP	●	●
SureTemp Plus Thermometry	●	●
Braun PRO 4000 Ear Thermometry	●	●
Masimo/Nellcor SpO ₂	●	●
Li-Ion Battery		
9-Cell (8 hours on-time performance)*	●	●
Communications		
USB	●	●
Ethernet	●	●
1 Host USB Accessory Port	●	○
4 Host USB Accessory Ports	○	●
Licenses		
Barcode	●	●
Spot Check Profile	●	●
Triage Profile	●	●
Weight Scale Connectivity**	●	●
Office Profile	●	○ (free)
Warranty		
2-Year Warranty	●	●

- Included
- Not Available
- Optional
- Upgradeable

Also available with
Braun ThermoScan[®] PRO
4000 Ear Thermometry



Connex VSM 6300—1 host USB port for accessories, 9-Cell Li-Ion Battery (featuring BP averaging)

63XXXX-B	CVSM 6300 – Blood Pressure
63XXPX-B	CVSM 6300 – Blood Pressure, Printer
63XTXX-B	CVSM 6300 – Blood Pressure, Temperature (SureTemp Plus)
63XTPX-B	CVSM 6300 – Blood Pressure, Temperature (SureTemp Plus), Printer
63XXE-B	CVSM 6300 – Blood Pressure, Temperature (Braun ThermoScan PRO 4000)
63XXPE-B	CVSM 6300 – Blood Pressure, Temperature (Braun ThermoScan PRO 4000), Printer
63NXXX-B	CVSM 6300 – Blood Pressure, SpO ₂ (Nellcor)
63NXPX-B	CVSM 6300 – Blood Pressure, SpO ₂ (Nellcor), Printer
63MXXX-B	CVSM 6300 – Blood Pressure, SpO ₂ (Masimo)
63MXPX-B	CVSM 6300 – Blood Pressure, SpO ₂ (Masimo), Printer
63NTXX-B	CVSM 6300 – Blood Pressure, SpO ₂ (Nellcor), Temperature (SureTemp Plus)
63NTPX-B	CVSM 6300 – Blood Pressure, SpO ₂ (Nellcor), Temperature (SureTemp Plus), Printer
63NXE-B	CVSM 6300 – Blood Pressure, SpO ₂ (Nellcor), Temperature (Braun ThermoScan PRO 4000)
63NXPE-B	CVSM 6300 – Blood Pressure, SpO ₂ (Nellcor), Temperature (Braun ThermoScan PRO 4000), Printer
63MTXX-B	CVSM 6300 – Blood Pressure, SpO ₂ (Masimo), Temperature (SureTemp Plus)
63MTPX-B	CVSM 6300 – Blood Pressure, SpO ₂ (Masimo), Temperature (SureTemp Plus), Printer
63MXXE-B	CVSM 6300 – Blood Pressure, SpO ₂ (Masimo), Temperature (Braun ThermoScan PRO 4000)
63MXPE-B	CVSM 6300 – Blood Pressure, SpO ₂ (Masimo), Temperature (Braun ThermoScan PRO 4000), Printer

Connex VSM 6400—4 host USB ports for accessories, 9-Cell Li-Ion Battery

64XXXX-B	CVSM 6400 – Blood Pressure
64XXPX-B	CVSM 6400 – Blood Pressure, Printer
64XTXX-B	CVSM 6400 – Blood Pressure, Temperature (SureTemp Plus)
64XTPX-B	CVSM 6400 – Blood Pressure, Temperature (SureTemp Plus), Printer
64XXE-B	CVSM 6400 – Blood Pressure, Temperature (Braun ThermoScan PRO 4000)
64XXPE-B	CVSM 6400 – Blood Pressure, Temperature (Braun ThermoScan PRO 4000), Printer
64NXXX-B	CVSM 6400 – Blood Pressure, SpO ₂ (Nellcor)
64NXPX-B	CVSM 6400 – Blood Pressure, SpO ₂ (Nellcor), Printer
64MXXX-B	CVSM 6400 – Blood Pressure, SpO ₂ (Masimo)
64MXPX-B	CVSM 6400 – Blood Pressure, SpO ₂ (Masimo), Printer
64NTXX-B	CVSM 6400 – Blood Pressure, SpO ₂ (Nellcor), Temperature (SureTemp Plus)
64NTPX-B	CVSM 6400 – Blood Pressure, SpO ₂ (Nellcor), Temperature (SureTemp Plus), Printer
64NXE-B	CVSM 6400 – Blood Pressure, SpO ₂ (Nellcor), Temperature (Braun ThermoScan PRO 4000)
64NXPE-B	CVSM 6400 – Blood Pressure, SpO ₂ (Nellcor), Temperature (Braun ThermoScan PRO 4000), Printer
64MTXX-B	CVSM 6400 – Blood Pressure, SpO ₂ (Masimo), Temperature (SureTemp Plus)
64MTPX-B	CVSM 6400 – Blood Pressure, SpO ₂ (Masimo), Temperature (SureTemp Plus), Printer
64MXXE-B	CVSM 6400 – Blood Pressure, SpO ₂ (Masimo), Temperature (Braun ThermoScan PRO 4000)
64MXPE-B	CVSM 6400 – Blood Pressure, SpO ₂ (Masimo), Temperature (Braun ThermoScan PRO 4000), Printer

Connex VSM Accessories

4800-60	Connex Accessory Cable Management Stand	36000	Braun PRO 4000 Dock (without Thermometer)
4701-62	Wall Mount with Basket	86100	Braun PRO 4000 Dock Wall-Mount Kit
6000-100H	Carrying Case, Hard Sided, VSM 6000	86200	Braun PRO 4000 Economy Mobile Stand-Mount Kit
6000-100S	Carrying Case, Soft Sided, VSM 6000	S1-6000	Connex VSM 6000 Comprehensive Partnership Program, 1-year agreement
103372	Model 6300 Spot Check Profile Authorization Code	S1-6000-2	Connex VSM 6000 Comprehensive Partnership Program, 2-year agreement
103373	Model 6300 Triage Profile Authorization Code	S2-6000	Connex VSM 6000 Biomed Partnership Program, 1-year agreement
4500-925	USB Cable for Wired Connectivity	S2-6000-2	Connex VSM 6000 Biomed Partnership Program, 2-year agreement
04000-900	Braun PRO 4000 Dock with Thermometer		

For a Connex VSM Office Profile license, please contact customer service at 1.800.535.6663

*On-time performance is an estimate based on using a new battery in a typical clinical use environment. Actual performance can vary based on workflow and configuration.

**Please contact your Health o meter representative for specific scale information.




Only Masimo® sensors should be used with Masimo Oximetry and Co-Oximetry. Masimo is a registered trademark of Masimo Corporation.



Only Nellcor® pulse oximetry sensors, including *OxiMax*® patient-dedicated adhesive sensors, should be used with the Nellcor pulse oximetry option. Nellcor and *OxiMax* are registered trademarks of Nellcor Puritan Bennett, Inc.



Braun® and ThermoScan® are registered trademarks of Braun, GmbH.



For more information about the Connex Vital Signs Monitor, please contact your local Welch Allyn representative or call 1.800.535.6663 to speak with a customer service representative.

Welch Allyn Corporate Headquarters
4341 State Street Road, P.O. Box 220, Skaneateles Falls, NY 13153-0220 USA
(p) 800.535.6663 (f) 315.685.3361 | www.welchallyn.com

©2012 Welch Allyn SM4088 Rev A

WelchAllyn®

Advancing Frontline Care™