

Cutimed® Sorbion®

Closing wounds. Together.

Cutimed® Sorbion®

Super-absorbent
wound dressing

For the
treatment of
highly exuding
wounds



THERAPIES. HAND IN HAND.

BSN medical

CUTIMED® SORBION®

The super-absorbent wound dressings of **Cutimed® Sorbion®** provide powerful exudate handling capacity. At the center of the unique design concept is patented Hydration Response Technology that comprises of a blend of harmonised product components – the unique fibre matrix, the gel forming polymers and the ultrasonically sealed casing ensure a secure treatment for highly exuding wounds.




Cutimed® Sorbion® provides powerful exudate management and outstanding retention capacity.



Fewer dressing changes = lower treatment costs




Convincing value



Optimal moist wound environment of the entire wound bed

- Ideal balance between exudate and maintaining a moist wound environment
- Gel forming polymers expand towards the wound, allowing for maximum wound contact even when highly saturated



Ideal wound bed preparation

- Soft debridement: The proportion of fibrinous coverings is reduced
- Binds inflammatory substances safely in the core
- Vertical absorption: protection against maceration of the wound edge



High wearing comfort

- Expansion edge technology allows **Cutimed® Sorbion®** to expand flat and evenly making it ideal to be used under compression
- Fibre matrix retains its shape even under saturation without the formation of lumps
- Dressing remains soft, flexible and conformable thanks to gel forming polymers
- Reduction of odours for better patient satisfaction


Performance at all levels

Indications

The **Cutimed® Sorbion®** range is indicated for use in highly exuding wounds such as:


- Pressure injuries
- Diabetic foot ulcers
- Venous leg ulcers

PERFORMANCE WITHOUT COMPROMISE



A DIVERSE RANGE OF ANATOMICALLY SHAPED DRESSINGS

Cutimed® Sorbion® offers a total range of diverse shapes and sizes to suit any wound, even on highly contoured parts of the body.



Cutimed® Sorbion® Sachet Multi Star
Cutimed® Sorbion® Sana Multi Star



Cutimed® Sorbion® Sachet S Drainage



Cutimed® Sorbion® Sachet XL

ADVANCED HEALING WITH CUTIMED® SORBION®



Ref. No.	Size	Pad Size	Items per Unit
Cutimed® Sorbion®			
7323200	7.5cm x 7.5cm		10
7323206	10cm x 10cm		10
7323209	20cm x 10cm		10
7323215	20cm x 20cm		10
7323218	30cm x 20cm		10
Cutimed® Sorbion® Sachet Multi Star			
7323701	ø 8cm		10
7323703	ø 14cm		10
Cutimed® Sorbion® Sachet XL			
7324000	45cm x 25cm		10
Cutimed® Sorbion® Sachet Border			
7323600	10cm x 10cm	7.5cm x 7.5cm	10
7323602	15cm x 15cm	10cm x 10cm	10
7323604	25cm x 15cm	20cm x 10cm	10
Cutimed® Sorbion® Sachet S Drainage			
7324100	10cm x 10cm		10
Cutimed® Sorbion® Sana Gentle			
7323319	8.5cm x 8.5cm		10
7323322	12cm x 12cm		10
7323325	22cm x 12cm		10
7323314	22cm x 22cm		10
7323316	32cm x 22cm		10
Cutimed® Sorbion® Sana Multi Star			
7323801	ø 11cm		10
7323803	ø 17cm		10

BSN medical Australia & New Zealand

L6, 211 Wellington Road,
Mulgrave, Victoria, Australia 3170

www.bsnmedical.com.au
www.bsnmedical.co.nz

BSN0145

10/2016

Customer Service Australia
T 1300 276 633 F 1300 998 830

Customer Service New Zealand
T 0508 276 111 F 0508 998 830

®Registered Trademark

