



Hillrom™

# Welch Allyn® Diagnostic Cardiology Suite™ Spirometry

## TECHNICAL SPECIFICATIONS

<b>Test Types Supported</b>	FVC, SCV, MVV, Pre- and Post-bronchodilation
<b>Sensor Type</b>	Pneumotach
<b>Accuracy &amp; Reproducibility</b>	Compliant with ATS/ERS 2005 Guidelines and ISO 26782 Guidelines
<b>Flow Range</b>	+/- 14 L/sec
<b>Volume Range</b>	0-14 L
<b>Resistance</b>	Compliant with ATS/ERS 2005 Guideline and ISO 26782 Guidelines (Impedance does not exceed 0.15kPa/(L/s) with flows up to 14 L/s)
<b>Interpretation</b>	1991 ATS Interpretation Standards Lung Age calculation available Automatic Interpretation available User-defined interpretation statements available
<b>Parameters</b>	FVC, FIVC, FIV1, FIV1%, FEV0.5, FEV1, FEV2, FEV3, FEV5, FEV6, FEV0.5, FEV0.5%, FEV1%, FEV1/FVC, FEV2%, FEV3%, FEV5%, FEV6%, PEF, FEF25, FEF50, FEF75, FEF0.2-1.2, FEF25-75, FEF75-85, PIF, FIF50, FEF50/FIF50, FEV1/FEV6, FET, MEF25, MEF50, MEF75 SVC, ERV, IRV, VT, IC, BF, MV, Tin, Tex, Tin/Tex MVV, MV, VT, BF, DFRC
<b>Quality Checks</b>	ATS Acceptability and ATS Reproducibility Effort-quality messages Audio and visual incentive for coaching patients

## MINIMUM SERVER REQUIREMENTS

<b>Processor</b>	Performance equivalent to an Intel Xeon® class, quad-core with hyperthreading
<b>Resolution</b>	1024×768 or better
<b>RAM</b>	Minimum 4 GB
<b>Operating System</b>	Microsoft Windows 2012 R2 (64-bit) Microsoft Windows 2016
<b>System Disk</b>	100 GB for OS and product installation (RAID recommended for data redundancy)
<b>Data Disks</b>	550 GB hard drive space available HD controller with 128 MB read/write cache (RAID recommended for data redundancy)



## MINIMUM STANDALONE/CLIENT SPECIFICATIONS

<b>Operating Systems</b>	Microsoft® Windows® 7 Professional (32-bit* or 64-bit) Microsoft Windows 10 Professional (64-bit) (Note: If installing ELI® Link with Diagnostic Cardiology Suite, system requires 64-bit OS to operate both.)	
<b>Compatible Viewing Resolutions and Font Scaling</b>	<b>Resolutions</b>	<b>Font Scaling</b>
	1024×768, 1366×768, 1280×800	100%
	1680×1050, 1920×1080, 1920×1200	100%, 125%, 150%
	2736×1824	100%, 125%, 150%, 200%
<b>Processor</b>	Performance equivalent or better than an Intel® Core™ i3	
<b>RAM</b>	Minimum 4 GB	
<b>Hard Drive Capacity</b>	Minimum 250 GB (If client: minimum 160 GB)	
<b>Archive</b>	Network or external USB drive	
<b>Software Installation</b>	USB	
<b>Network</b>	100 Mbps connection or better	
<b>USB Ports</b>	Requires at least two USB 2.0 ports (spirometer external USB drive)	
<b>Input Devices</b>	Standard keyboard, mouse, spirometer	

## MINIMUM CITRIX® XENAPP® REQUIREMENTS

<b>Client Machines That Run Citrix XenApp</b>	Microsoft Windows 7 Professional (64-bit and 32-bit) Microsoft Windows 10 Professional (64-bit) Citrix Receiver Any Internet browser supported by Citrix
<b>Citrix Desktop Delivery Controller</b>	Citrix XenApp version 7.13 Any operating system supported by Citrix
<b>Citrix App Server</b>	Microsoft Windows 2012 R2 (64-bit) Microsoft Windows 2016 Citrix Virtual Delivery Agent 7.13

## MINIMUM WINDOWS RDP REQUIREMENTS

<b>Thin Clients Using Windows RDP</b>	Microsoft Windows 7 Professional (64-bit and 32-bit) Microsoft Windows 10 Professional (64-bit)
<b>Operating System for Servers Running Windows RDP</b>	Microsoft Windows 2012 R2 (64-bit) Microsoft Windows 2016

Contact your Hillrom representative for more information.



hillrom.com

\*If using EMR interface, font scaling may vary.  
Specifications are subject to change without notice.