

Rainbow Trays™

Colour coded Anaesthetic trays





- Hypnotics**
 - Diazepam**mg/ml.
 - Midazolam**mg/ml.
- Opioids**
 - Morphine**mg/ml.
 - Fentanyl**micrograms/ml.
- Induction agents**
 - Propofol**mg/ml.
 - Ketamine**mg/ml.
- Neuromuscular blocking drugs**
 - Suxamethonium**mg/ml.
 - Vecuronium**mg/ml.
- Anti-emetics**
 - Ondansetron**mg/ml.
- Miscellaneous**
 - Heparin**units/ml.
 - Protamine**mg/ml.
- Neuromuscular blocking drug antagonist**
 - Neostigmine**micrograms/ml.
- Vasopressors**
 - Adrenaline**micrograms/ml.
 - Ephedrine**mg/ml.
- Anticholinergic agents**
 - Atropine**micrograms/ml.
 - Glycopyrronium**micrograms/ml.
- Local anaesthetics**
 - Lidocaine**%.
 - Bupivacaine**%.

Rainbow Trays™

Colour coded Anaesthetic trays




RAINBOW TRAYS™ BASES & INSERTS




Inserts locate in one orientation only

LOADING THE TRAYS




Syringe labels are colour coded to each section

RECYCLING THE TRAYS



The base trays are reusable and can be anti-bacterially wiped clean or sterilised with Ethylene Oxide



The single use insert trays are 100% recyclable