

NIO-Pediatric

Quick, Safe, Easy Vascular Access

- Vascular access in under 10 seconds
- Adjustable according to patient age
- Redundant safety features
- Semi-automatic device
- No exposed needle
- Single use device
- 5 year shelf life
- Self-contained



Vascular Access in a Flash

When you need immediate vascular access, you can be confident that the NIO-P is complete and ready to use. The NIO-P is the quickest, most reliable intraosseous device for gaining vascular access in pediatric patients.

Product Highlights

Semi-automatic, self-contained unit facilitates rapid vascular access in under 10 seconds.

Unexposed needle decreases the risk of accidental needle injury.

Redundant safety features prevent accidental activation of the device.



Specifications

Dimensions	148.55 mm (5.9") L 64.5 mm (2.6") W 43.7 mm (1.72") H
Weight	85.6 g
Packaging	Packed individually in a peel back blister pack
Needle Gauge	18 G
Needle Effective Length	38.1 mm (1.5")
Needle Penetration Depth	(Age 9-12): Closed stabilizer - 18 mm (Age 3-9): Open stabilizer - 14 mm
Safety Mechanism	1. Rotating safety latch prevents accidental activation 2. Simultaneous pressing of device against the bone and pulling of trigger mechanism required to activate device
Shelf Life*	5 years

Specifications are subject to change without notice.
*From date of manufacture

The NIO-P features an innovative stabilizer hub which serves four functions:

- **Unique location arrows** aid insertion site identification
- **Needle stabilizer** assists in securing the needle to the limb.
- **Adjustable needle penetration depth** means one device can accommodate a wide range of pediatric patients
- **Mechanical stopper** limits needle depth, preventing over penetration

PerSys Medical
5310 Elm St.
Houston, TX 77081 USA
1.888.737.7978
info@ps-med.com
www.ps-med.com

