

DONJOY®

DonJoy
StrapiLax™ Knee



DonJoy StrapiLax™ Knee

Bringing together flexible bilateral stays, a contoured knee pad, and adjustable strapping, DonJoy StrapiLax™ Knee helps relieve, stabilize, and protect the knee joint from abnormal movements, as well as aiding proprioceptive acuity.

Intended use/Indications

- Support and compression following acute and chronic knee injury and mild knee osteoarthritis
- Proprioceptive and neuromuscular stimulation through compression
- Mild medial (inner) and lateral (outer) ligament sprains of the knee
- Mild knee instability
- Prevention of injury

Characteristics

- Knitted construction using a multidirectional elastic fabric
- Lightweight
- Low profile
- Stylish, modern design
- Fits left or right knee
- Not made with natural rubber latex
- Easy to clean, machine washable
- Available in pediatric size

Features and benefits



Anatomical 3D knit design

to help provide effective and targeted compression (20-36mmHg, medical standard) and support



Breathable, elastic, soft knit

for enhanced comfort



Stretch zone

over the back of the knee to aid breathability, comfort, and freedom of movement



Compression-reduced edges

help divert pressure at the ends of the support to provide a comfortable fit



Anatomically contoured knee pad

to aid stability of the patella, massaging effect, and to help improve proprioception



Non-slip silicone-coated bands

help provide extra secure hold under stress



Adjustable knee strapping

works like functional tape to help provide added compression, stability, and a custom fit



ORDERING INFORMATION

PART NUMBER	KNEE CIRCUMFERENCE		SIZE
	CM	Inches	
82-0711-0	30 - 34 cm	11¾ - 13½"	0/PEDIATRIC
82-0711-1	35 - 38 cm	13¾ - 15"	1
82-0711-2	39 - 42 cm	15½ - 16½"	2

82-0711-3	43 - 46 cm	17 - 18"	3
82-0711-4	47 - 49 cm	18½ - 19¼"	4
82-0711-5	50 - 53 cm	19¾ - 20 ¾"	5
82-0711-6	54 - 57 cm	21¼ - 22½"	6
82-0711-7	58 - 61 cm	22¾ - 24"	7

Measurement taken 15 cm above mid-patella



DJO Global
1a Guildford Business Park
Guildford, Surrey
GU2 8XG
UNITED KINGDOM

djoglobal.eu

© 2021 DJO - 000-00001372-10-EN - Rev A
Individual results may vary. Neither DJO Global, Inc. nor any of its subsidiaries dispense medical advice. The contents of this brochure do not constitute medical, legal, or any other type of professional advice. Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you.