

BARRIER® Surgical gown

Primary Plus





Shoulder pleats allow freedom of movement for the upper body




Ergonomic sleeves to allow comfort and freedom of movement

- Extra space in shoulder to provide freedom of movement
- Hook and loop secured to the nonwoven
- Meets EN 13795-1:2019
- Meet AAMI 2 (SP)
- Meet AAMI 3 (HP)


Standard performance	Sterile 2 towels and wrapped	Sterile	Non-sterile (tray code)	Size	Length	Width	Disp/Trp
	580200-00	580100-00	480100-00 (2325283-00)	L	119 cm	150 cm	20/40
	580201-00	580101-00	480101-00 (2325269-00)	L-L	134 cm	150 cm	18/36
	580202-00	580102-00	480102-00 (2325246-00)	XL	129 cm	170 cm	18/36
	580203-00	580103-00	480103-00 (2325290-00)	XL-L	144 cm	170 cm	16/32
	580204-00	580104-00	480104-00 (2325291-00)	XXL-L	145 cm	184 cm	16/32
	580205-00	580105-00	480105-00 (2325284-00)	XXL-XL	157 cm	184 cm	16/32

High performance	Sterile 2 towels and wrapped	Sterile	Non-sterile (tray code)	Size	Length	Width	Disp/Trp
	590200-00	590100-00	490100-00 (2325247-00)	L	119 cm	150 cm	18/36
	590201-00	590101-00	490101-00 (2325274-00)	L-L	134 cm	150 cm	16/32
	590202-00	590102-00	490102-00 (2325248-00)	XL	129 cm	170 cm	16/32
	590203-00	590103-00	490103-00 (2325292-00)	XL-L	144 cm	170 cm	14/28
	590204-00	590104-00	490104-00 (2325237-00)	XXL-L	145 cm	184 cm	14/28
	590205-00	590105-00	490105-00 (2325293-00)	XXL-XL	157 cm	184 cm	14/28



- Impermeable sleeve reinforcement to minimise risk of strike through
- Impermeable front reinforcement to minimise risk of strike through
- Conformable material
- Front reinforcement is breathable

Length: 119cm



Circumference: 150cm

Measurements for size L

Product classification

Standard performance gowns: European Standard EN 13795 Standard performance and ANSI/AAMI PB70 level 2. High performance gowns: European Standard EN 13795 High Performance and ANSI/AAMI PB70 level 3

Packaging

Three layers of packaging, sterility barrier, retail box and transport box, help reduce risk of contamination and damage during transportation and storage.