



# USL STABILISER MOON BOOT

Maximum stability and support following injury  
or post-operative procedures

## CONTOURED STRUT DESIGN

Allows the walker  
to conform to the leg,  
ensuring a better fit.

**PLASTIC MOULDED  
UPRIGHTS** with steel  
reinforcement, for increased  
durability.

## ROCKER SOLE

Helps promote a natural  
gait to allow continuation  
of daily activities.



## WIDE FOOTBED

Provides increased  
stability, comfort and  
accommodation for edema  
and bandages.

## SHOCK- ABSORBING INSOLE

Reduces impact of  
heel strike when  
walking.



CODE	DESCRIPTION
------	-------------

22149	USL Stabiliser Tall Moon Boot X-Small
22150	USL Stabiliser Tall Moon Boot Small
22151	USL Stabiliser Tall Moon Boot Medium
22152	USL Stabiliser Tall Moon Boot Large
22153	USL Stabiliser Tall Moon Boot X-Large


# Indications

- Acute ankle sprains
- Soft tissue injuries
- Stress fractures of the lower leg
- Stable fractures of the foot and ankle
- Achilles tendon repairs
- Post-operative procedures



# Find your size

BOOT SIZE	MENS SHOE SIZE	WOMENS SHOE SIZE	FOOT LENGTH
X-Small	Kid's 1½ - Men's 4	Kid's 1½ - Women's 5½	20 - 22cm
Small	4½ - 6½	6 - 8	22 - 24cm
Medium	7 - 9	8½ - 10½	24.5 - 26cm
Large	9½ - 12	11 - 12+	26.5 - 28.5cm
X-Large	12½ - 14		29 - 30cm



For foot length,  
measure foot from  
heel to toe.

## SPECIAL CONSIDERATIONS



Consider a **LARGER BOOT** if you're near the maximum size **AND** these apply:

- You wear wide shoes
- You will use bandages or thick socks
- You have swelling
- You typically like a looser fit



Consider a **SMALLER BOOT** if you're near the minimum size **AND** these apply:

- You wear narrow shoes or have slim feet
- You typically like a tighter fit



**JO RICE**  
Product Manager -  
Braces & Support

visit our website  
**uslmedical.co.nz**

order by freephone  
**0800 658 814**

customerservices@uslmedical.co.nz

