

SAFETY DATA SHEET

Section 1.	Identification of the material and the supplier
-------------------	--

Product: **USL Lubrigel**
 Item Code:
 Product Use: Lubricant (Non Sterile)
 Restriction of Use: External use only

New Zealand Manufacturer: **Jaychem Industries Ltd**
 Address: **3 Kordel Place
 East Tamaki, Auckland**

Telephone: +64 9 274 6647
 Fax Number: +64 9 274 1358

Emergency Telephone: 0800 764 766 (National Poison Centre)

Date of SDS Preparation: 1 April 2021

Section 2.	Hazards Identification
-------------------	-------------------------------

This substance is hazardous according to the HSNO (Minimum Degrees of Hazard) Regulations 2001

EPA Approval No: Cosmetic Group Standard – HSR002552

Pictograms



Signal Word: Warning

HSNO Classification	Hazard Code	Hazard Statement	GHS Category
6.5B	H317	May cause an allergic skin reaction.	Category 1

Prevention Code	Prevention Statement
P103	Read label before use.
P261	Avoid breathing vapours.
P272	Contaminated work clothing should not be allowed out of the workplace.
P280	Wear protective clothing.

Response Code	Response Statement
P363	Wash contaminated clothing before reuse.
P302 + P352	IF ON SKIN: Wash with plenty of soap and water.
P333 + P313	If skin irritation or rash occurs: Get medical advice/attention.

Storage Code	Storage Statement

None allocated	
----------------	--

Disposal Code	Disposal Statement
P501	Dispose of according to Local Regulations or Authorities

Section 3. Composition / Information on Ingredients

Ingredients	Wt%	CAS NUMBER.
Hydroxy propylmethyl cellulose	<2	9004-65-3
Preservative	0.5-1.5	Proprietary
Non-hazardous ingredients and Deionised Water	Balance	

Section 4. First Aid Measures

Routes of Exposure:

- If in Eyes Rinse cautiously with water for 15 minutes. If eye irritation persists: Get medical advice.
- If on Skin Wash with plenty of soap and water. If skin irritation or rash occurs: get medical advice/attention.
- If Swallowed Rinse mouth. Never give anything to the mouth of an unconscious person. If vomiting occurs, place victim face downwards, with the head turned to the side and lower than the hips to prevent vomit entering the lungs. Seek medical attention if needed.
- If Inhaled Remove person to fresh air. Remove contaminated clothing and loosen remaining clothing. Allow person to assume most comfortable position and keep warm. Keep at rest until fully recovered. Get medical advice if breathing becomes difficult.

Section 5. Fire Fighting Measures

Hazard Type	Non Flammable.
Hazards from combustion products	None known.
Suitable Extinguishing media	Use media suitable for surrounding materials.
Precautions for firefighters and special protective clothing	Wear protective gear.
HAZCHEM CODE	None Allocated

Section 6. Accidental Release Measures

- For bulk spillages:**
Wear protective equipment as detailed in Section 8.
- Scrape up and place in a container for disposal.
- Dispose of according to Local Regulations.

Section 7. Handling and Storage

Precautions for Handling in bulk:

Product Name: USL Lubrigel
Date of SDS: 1 April 2021

Issued by: Technical Compliance Consultants (NZ) Ltd
Tel: 64 9 475 5240 www.techcomp.co.nz

- Read label before use.
- Avoid breathing vapours.
- Contaminated work clothing should not be allowed out of the workplace.
- Wear protective clothing.

Precautions for Storage in bulk:

- Store away from incompatible materials listed in Section 10.
- Keep containers securely sealed.
- Store in a cool, dry and well-ventilated area.
- Protect containers against physical damage.
- Check regularly for spills and leaks.

Section 8 Exposure Controls / Personal Protection

WORKPLACE EXPOSURE STANDARDS (provided for guidance only)

Substance	TWA		STEL	
	ppm	mg/m ³	ppm	mg/m ³
Propane-1,2 diol [57-55-6] vapour & particulates	150	474	-	-

Workplace Exposure Standard – Time Weighted Average (WES-TWA). *The time-weighted average exposure standard designed to protect the worker from the effects of long-term exposure.* Workplace Exposure Standard – Short-Term Exposure Limit (WESSTEL). *The 15-minute average exposure standard.* Applies to any 15- Minute period in the working day and is designed to protect the worker against adverse effects of irritation, chronic or irreversible tissue change, or narcosis that may increase the likelihood of accidents. The WES-STEL is not an alternative to the WES-TWA; both the short-term and time-weighted average exposures apply.

Engineering Controls

Ensure that ventilation maintains dust levels below WES.

Personal Protection Equipment

Eyes	Not required for normal use.
Hands and Skin	Not required for normal use.
Respiratory	Not required for normal use.

Section 9 Physical and Chemical Properties

Appearance	Clear colourless viscous liquid
Odour	Not available
Odour Threshold	Not applicable
pH	5.6
Boiling Point	Not applicable
Melting Point	Not applicable
Freezing Point	Not applicable
Flash Point	Not applicable
Flammability	Not applicable
Upper and Lower Explosive Limits	Not applicable
Vapour Pressure	Not applicable
Vapour Density	Not applicable
Specific Gravity	0.95 – 1.05
Solubilities	Not applicable
Partition Coefficient:	Not applicable
Auto-ignition Temperature	Not applicable
Decomposition Temperature	Not applicable

Kinematic Viscosity	Not applicable
Particle Characteristics	Not applicable

Section 10. Stability and Reactivity

Stability of Substance	This product is stable under normal conditions.
Conditions to Avoid	Avoid freezing.
Incompatible Materials	None known.
Hazardous Decomposition Products	No hazardous products are expected.

Section 11 Toxicological Information

Acute Effects:

Swallowed	Not applicable.
Dermal	Not applicable.
Inhalation	Not applicable.
Eye	Not applicable.
Skin	May cause an allergic skin reaction.

Chronic Effects:

Carcinogenicity	Not applicable.
Reproductive Toxicity	Not applicable.
Germ Cell Mutagenicity	Not applicable.
Aspiration	Not applicable.
STOT/SE	Not applicable.
STOT/RE	Not applicable.

Section 12. Ecotoxicological Information

This product is not known to be a hazard to the environment.

Persistence and degradability	No data available
Bioaccumulation	No data available
Mobility in Soil	No data available
Other adverse effects	No data available

Section 13. Disposal Considerations

Disposal Method: Dispose of product and packaging at an approved landfill or other approved Facility.

Precautions: Preferably dispose of by use.

Disposal methods to avoid: None known.

Section 14 Transport Information

This product is not classified as a Dangerous Good for transport in NZ ; NZS 5433:2012

Section 15 Regulatory Information

EPA Approval Code: Cosmetic Products Group Standard – HSR002552

HSNO Classification: 6.5B

Product Name: USL Lubrigel
Date of SDS: 1 April 2021

Issued by: Technical Compliance Consultants (NZ) Ltd
Tel: 64 9 475 5240 www.techcomp.co.nz

HSNO Controls:
Trigger quantities for this substance:

	Trigger Quantity
Approved Handler	Not required
Location Certificate	Not required
Tracking Trigger Quantities	Not required
Signage Trigger Quantities	Not required
Emergency Response Plan	Not required
Secondary Containment	Not required
Restriction of Use	For external use only.

Section 16 Other Information

Glossary

EC ₅₀	Median effective concentration.
EEL	Environmental Exposure Limit.
EPA	Environmental Protection Authority
HSNO	Hazardous Substances and New Organisms.
LC ₅₀	Lethal concentration that will kill 50% of the test organisms inhaling or ingesting it.
LD ₅₀	Lethal dose to kill 50% of test animals/organisms.
LEL	Lower explosive level.
OSHA	American Occupational Safety and Health Administration.
TEL	Tolerable Exposure Limit.
TLV	Threshold Limit Value-an exposure limit set by responsible authority.
UEL	Upper Explosive Level
WES	Workplace Exposure Limit

1. HSNO Approved Code of Practice: Preparation of Safety Data Sheets, September 2006.

Disclaimer

This document has been issued by TCC (NZ) Ltd and serves as their Safety Data Sheet ('SDS'). It is based on information concerning the product which has been provided to TCC (NZ) Ltd or obtained from third party sources and is believed to represent the current state of knowledge as to the appropriate safety and handling precautions for the product at the time of issue. Further clarification regarding any aspect of the product should be obtained directly from the manufacturer. While TCC (NZ) have taken all due care to include accurate and up-to-date information in this SDS, it does not provide any warranty as to accuracy or completeness. As far as lawfully possible, TCC (NZ) Ltd accept no liability for any loss, injury or damage (including consequential loss) which may be suffered or incurred by any person as a consequence of their reliance on the information contained in this SDS

The information herein is given in good faith, but no warranty, express or implied is made.

Please contact the New Zealand distributor, if further information is required.

Issue Date: 1 April 2021 Review Date: 1 April 2026