

Reynard Premier biodegradable detergent & disinfectant wipes

A single step solution to cleaning and disinfecting surfaces and equipment, proven to kill at least 99.999% of pathogens. With a superior cloth and unsurpassed equipment compatibility, germs don't stand a chance.

Specifications

- Effective for use against COVID and SARS
- Unparalleled equipment compatibility
- Superior cloth: large size versus most competing wipes

Points to remember

- Use one wipe per piece of equipment or surface
- Wipe from clean to dirty
- Clean from top to bottom

Testing Data.

All claims in accordance with TGO 104 Order 2019

Effective against a wide range of pathogenic micro-organisms including bacteria, viruses such as SARS-CoV-1&2 and MERS-CoV virus fungi and spores.

PATHOGEN

EN Test

Kills at least 99.999% of pathogens. Effective from 60 seconds

Bacteria	EN13727 and EN16615: 2015
Fungi	EN13624 and EN16615: 2015
Sporicidal	EN13704 and EN17126
Microplasma	EN1276
Viruses	EN14476
Corona Virus	EN14476

- **Biodegradable**
- **Made with 100% Bamboo**



SKU	WIPES PER PACK	PACKS PER CARTON
RHS210 - Biodegradable	50	12



Reynard is reassurance.



Australia

44 Hampden Road
Artarmon, NSW, 2064
enquiries@reynardhealth.com.au
1300 667 199

New Zealand

PO Box 8470
Havelock North
Sales@reynardhealth.com
0508 494 737

UK

13 White Horse Business Park,
Stanford in the vale SN7 8NY
uk@reynardhealth.com
+44 (0)845 263 7335

International

PO Box 8470
Havelock North
Sales@reynardhealth.com
+64 (0)6 650 0709