



Hillrom™



Welch Allyn
SureTemp Plus
Electronic Thermometers



WELCH ALLYN SURETEMP PLUS 690 ELECTRONIC THERMOMETER

- Removable probe well for easy cleaning
- Large, easy-to-read LCD
- C/F button for toggling between Fahrenheit and Celsius
- Recall button displays last temperature taken
- Rubberized grips
- Convenient storage for 25 probe covers
- Battery life indicator
- Mode button



OPTIONAL WALL MOUNT

The SureTemp Plus 690 with optional wall mount has storage for an additional 50 probe covers.

WELCH ALLYN SURETEMP PLUS 692 ELECTRONIC THERMOMETER

SAME FEATURES AS THE 690 MODEL, PLUS:

- NEW** Backlit display helps in low-light environments
- Pulse timer button/monitor mode

NEW FEATURE Easier to switch to Monitor Mode for measuring a patient's thermal steady state (approx. 3 to 5 minutes)



STANDARD WALL MOUNT

The SureTemp Plus 692 with standard wall mount has storage for an additional 50 probe covers.

Interchangeable rectal probe and well reduce risk of cross-contamination

IDEAL CARE SETTINGS:

- Family Practice
- Pediatrics
- Long-term Care
- Internal Medicine
- OB/Gyn



IDEAL CARE SETTINGS:

- Urgent Care
- MedSurge
- Maternity
- Emergency Medicine
- Oncology
- Cardiology

SETTING THE STANDARD OF CARE

Enhance a routine task with versatile, easy and accurate technology featuring oral, axillary and rectal modes to meet all your temperature taking needs.

1 EASY

An easy-to-read LCD display with flashing icons, and recall button displays last temperature taken. SureTemp Plus 692 model features a 60-second timer to support pulse and respiration readings.

3 FAST AND ACCURATE

The SureTemp Thermometer captures quick readings - approximately 4-6 seconds oral, 10-13 second pediatric axillary and rectal readings and 12-15 second adult axillary readings. Toggle between Celsius and Fahrenheit, or monitor mode that helps measure thermal steady state.

2 SAFE AND SECURE

Soft, single-use disposable covers help minimize cross-contamination risk, and interchangeable probe wells are easy to remove and clean as you switch between modes. The SureTemp Plus 692 model is designed with two security features to help protect against theft.



User-selectable icon and/or words for identification of thermometer mode

Battery life indicator

TECHNICAL SPECIFICATIONS

Calibration Accuracy	±0.2 °F (±0.1 °C) (Monitor Mode)
Temperature Predict Time	Oral: Approximately 4-6 seconds Adult Axillary (18 years and older): Approximately 12 - 15 seconds Pediatric Axillary (17 years and younger): Approximately 10 - 13 seconds Rectal: Approximately 10 - 13 seconds
Unit Dimensions	8.46" x 3.18" x 2.43" (215 mm x 81 mm x 62 mm)
Weight	Model 690: 12.6 oz (357 g), Model 692: 12.6 oz (357 g)
Ambient Temperature Range	50.0 °F (10 °C) to 104 °F (40.0 °C)
Patient Temperature Range	80.0 °F (26.7 °C) to 110 °F (43.3 °C)
Power Source	(3) 1.5 Vdc AA batteries
Battery Life	Model 690: Approximately 6,000 temperature, Model 692: Approximately 5,000 temperatures

PRODUCT AND ACCESSORY ORDERING INFORMATION

THERMOMETERS

690 MODELS:

01690-200: SureTemp Plus 690 Electronic Thermometer, Oral Probe and Oral Probe Well with 4ft cord
01690-201: SureTemp Plus 690 Electronic Thermometer, Rectal Probe with Rectal Probe Well with 4ft cord
01690-300: SureTemp Plus 690 Electronic Thermometer, One per room, Wall Mount, Oral Probe with Oral Probe Well with 9ft cord (must be locked to the wall)
01690-301: SureTemp Plus 690 Electronic Thermometer, One Per Room, Wall Mount, Rectal Probe with Rectal Probe Well with 9ft cord (must be locked to the wall)
01690-400: SureTemp Plus 690 Electronic Thermometer, Wall Mount, Oral Probe with Oral Probe Well with 4ft cord
01690-700: SureTemp Plus 690 Electronic Thermometer, Wall Mount, Oral Probe with Oral Probe Well with 9ft cord with Rolling Stand and Holder for Spare Probe & Well

ACCESSORIES

21326-0000: Wall Holder
21330-0000: Holder for Spare Probe and Well
06138-000: Calibration Key Assembly, 690/692
02892-000: Probe and Well Kit, 4ft, Rectal
02892-100: Probe and Well Kit, 9ft, Rectal
02893-000: Probe and Well Kit, 4ft, Oral
02893-100: Probe and Well Kit, 9ft, Oral
02891-0000: Probe Well, Blue

692 MODELS:

01692-200: SureTemp Plus 692 Electronic Thermometer, Wall Mount, Security System with ID Location Field, Oral Probe and Oral Probe Well with 4ft cord
01692-201: SureTemp Plus 692 Electronic Thermometer, Wall Mount, Security System with ID Location Field, Rectal Probe and Rectal Probe Well with 4ft cord
01692-300: SureTemp Plus 692 Electronic Thermometer, One per room, Wall Mount, Oral Probe and Oral Probe Well with 9ft cord (must be locked to the wall)
01692-301: SureTemp Plus 692 Electronic Thermometer, One per room, Wall Mount, Rectal Probe and Rectal Probe Well with 9ft cord (must be locked to the wall)
01692-700: SureTemp Plus 692 Electronic Thermometer, Wall Mount, Oral Probe and Oral Probe Well with 9ft cord, with Rolling Stand and Holder for Spare Probe & Well

02891-1000: Probe Well, Red
08273-000: Thermometer Stand with Holder for Spare Probe and Well
21333-0000: Bracket, M600 Stand Upgrade, with Holder for Spare Probe and Well
98025-0000: Extended Warranty M692

REPLACEMENT PROBE COVERS

05031-101: Case of 1,000 Probe Covers

For more information, please contact your local distributor or Hillrom sales representative at 1-800-535-6663.

hillrom.com

Baxter, Hillrom, SureTemp, and Welch Allyn are trademarks of Baxter International, Inc. or its subsidiaries.

US-FLC166-230001 V2 02/23