


	Circumference in cm	xx-small *	x-small	x-small PLUS	small	small PLUS	medium	medium PLUS	large	large PLUS
 cT	cT	60-75	65-80	75-90	75-90	85-100	85-100	95-110	95-110	105-120
cH	cH	75-90	80-95	85-100	85-100	95-110	95-110	105-120	105-120	115-130
cG	cG	42-54	47-59	52-64	52-64	56-68	56-68	61-74	61-74	66-79
cF	cG Grip top	42-52	47-56	52-61	52-61	56-65	56-65	61-71	61-71	66-76
cE	cF	38-49	41-52	44-55	44-55	49-60	49-60	53-64	53-64	59-70
cD	cE	29-36	33-40	37-44	37-44	39-46	39-46	41-48	41-48	45-52
cC	cD	23-30	27-34	31-38	31-38	34-41	34-41	37-44	37-44	41-48
cB1	cC	25-32	29-36	33-40	33-40	35-42	35-42	38-45	38-45	42-49
cB	cB1	22-26	24-29	24-29	26-32	26-32	29-35	29-35	31-37	31-37
cY	cB	18-20	20-22	20-22	22-24	22-24	24-26	24-26	26-29	26-29
cA	cY	27-31	28-33	28-33	29-35	29-35	31-37	31-37	32-38	32-38
	cA	18-20	20-22	20-22	22-24	22-24	24-26	24-26	26-29	26-29

*only in SKIN and BLACK, only satin knobbed grip top



Length in cm					
short	—	short	56–64	short	—
normal	35–40	normal	65–74	normal	65–74
long	41–45	long	75–84	long	75–84
					