



Zetuvit®

The ideal choice for managing low to moderately exuding wounds.

Zetuvit® is a multi-layered dressing pad for managing wound exudate. The high absorption and retention capacity makes Zetuvit® an economic choice, reducing the frequency for dressing changes, while assisting in protecting the wound and peri wound from maceration.

Zetuvit® is enclosed in a soft, non-woven hydrophobic polyamide fibre which helps prevent the dressing from adhering to the wound. The inner hydrophilic surface of the cover has a high capillary activity providing wicking of wound exudate into the next layer of cellulose tissue. The tissue layer rapidly disperses fluid into the central core of highly absorbent cellulose fluff where it is retained.

Zetuvit® has a backing layer of hydrophobic cellulose which features a special water repellent, non-woven material that is permeable to air and helps prevent strike through.

Zetuvit® can be used to manage low to moderately exuding wounds, as a primary dressing on surgical wounds or under compression bandaging. The frequency of dressing changes for Zetuvit® is dependent upon the amount of wound exudate.

Zetuvit® will require a retention tape such as Omnifix® or can be secured with a light cohesive bandage such as Peha-haft®.

Applications

- Manages low to moderate exuding wounds
- Assists with wound bed preparation by effective exudate management
- Suitable as a primary dressing on surgical wounds

Benefits

- Low irritant and low adherent
- Conformable
- Cushions wound for protection
- Protects the wound and peri wound from maceration
- Cost and time effective
- Water repellent top layer to prevent strike through



Product	Size	Unit	Code	Units Per case
Zetuvit®	10 x 10 cm	25 pcs	413 701	10
	10 x 20 cm	25 pcs	413 702	10
	20 x 20 cm	15 pcs	413 703	9
	20 x 40 cm	5 pcs	413 704	9
	13.5 x 25 cm	10 pcs	413 705	6

1800 805 839
www.hartmann.com.au



Zetuvit® Plus

The ideal choice for managing heavily exuding wounds.

Zetuvit® Plus is a multi-layered dressing pad for managing wound exudate and contains Super Absorbent Polymers (SAP) for rapid and reliable fluid retention.

Zetuvit® Plus has a four layer composition. The absorbent core features cellulose fluff bound to SAP, which rapidly binds wound exudate and reliably retains it within the absorbent core. The cellulose fluff also acts to cushion and protect the wound.

The absorbent core is enclosed in a thin, non-woven fabric that uniformly distributes the fluid, while the entire dressing is enclosed in a soft two-layer outer non-woven material. The soft hydrophobic outer surface of the non-woven polyamide helps prevent the dressing from adhering to the wound.

The green side of the dressing indicates the water repellent top layer which is a non-woven material that is permeable to air and helps prevent strike through

The high absorption and retention capacity makes Zetuvit® Plus an economic choice, reducing the frequency for dressing changes while assisting in protecting the wound and peri wound skin from maceration.

Zetuvit® Plus is an ideal dressing for difficult to manage, highly exuding wounds such as fungating tumours and wounds associated with lymphoedema.

Zetuvit® Plus can be used to manage moderate to severely exuding wounds. The frequency of dressing changes for Zetuvit® Plus is dependent upon the amount of wound exudate.

Zetuvit® Plus will require a retention tape such as Omnifix® or can be secured with a light cohesive bandage such as Peha-haft®.

Applications

- Manages moderate to heavy exuding wounds
- Assists with wound bed preparation by effective exudate management
- Suitable for use under compression bandages

Benefits

- Absorbent core contains SAP for rapid and reliable fluid retention
- Green water repellent top layer to prevent strike through
- Low irritant and low adherent
- Protects the wound and peri wound from maceration
- Cost and time effective



Product	Size	Unit	Code	Units Per case
Zetuvit® Plus	10 x 10 cm	10 pcs	413 710	8
	10 x 20 cm	10 pcs	413 711	6
	15 x 20 cm	10 pcs	413 712	8
	20 x 25 cm	10 pcs	413 713	6
	20 x 40 cm	5 pcs	413 714	8

1800 805 839
www.hartmann.com.au