

Omnifix[®]

The self-adhesive, polyacrylate non-woven fabric retention dressing

Omnifix elastic is a soft, self-adhesive fabric for secure retention of wound dressings. The hypoallergenic acrylic adhesive reduces the risk of allergic reaction.

Omnifix elastic is permeable to both air and water vapour, enabling large areas to be covered without the risk of maceration.

The widthways elasticity of Omnifix elastic allows its smooth application to joints and angular parts of the body, with no inhibition of movement.

Omnifix elastic covering paper has a 'wave' cut to aid in its removal and a measuring grid printed on it to facilitate cutting to size.

Omnifix E is also available as a cost-effective part of this range. This product has a non-woven polyester backing and a synthetic rubber based adhesive.





Product	Width	Unit	Code	Units Per case
Omnifix® elastic Individually boxed 2 m length 10 m length	10 cm	1 roll	900 601	24
	2.5 cm	2 rolls	900 606	60
	5 cm	1 roll	900 602	60
	10 cm	1 roll	900 603	32
	15 cm	1 roll	900 604	20
Omnifix® E In cardboard boxes 10 m length	20 cm	1 roll	900 605	16
	30 cm	1 roll	900 607	10
	5 cm	1 roll	900 649	60
	10 cm	1 roll	900 650	32
	15 cm	1 roll	900 651	20
	20 cm	1 roll	900 652	16