

RelevoTM

Relevo has the following product usage data:

Absorption capacity: Average value for free absorption of saline solution 0,9%:
1 ml per cm² of core size.

RELEVO

REF	Product Size	Free swell abs cap, saline solution 0,9 %	Packing: Pieces/ inner box	Packing: Inner boxes per outer box	CE class
F50011/50	10 x 10 cm	50 ml	50	9	IIb
F50011/10			10	45	IIb
STERILE				(450 pcs)	
F50021/50	10 x 20 cm	120 ml	50	8	IIb
F50021/10			10	45	lib
STERILE				(450 pcs)	
F50041/50	20 x 20 cm	240 ml	50	6	IIb
F50041/10			10	30	IIb
STERILE				(300 pcs)	

Use

Indications: Exuding wounds such as leg ulcers, pressure ulcers, diabetic foot, etc.

Special features: Absorbs more exudate than cellulose dressings, and
absorbs under compression.
Retains exudate, will not leak or drip during or after absorption.

QM/20140128