

# EZ100m



Single Use Atomizer and Power-sprayer. EZ100m was designed for the rapid anesthetization and/or spraying of patients using aqueous or viscous medications or solutions.

## EZ100m Features and Benefits

- 15 cc calibrated clear bottle allows for metered dose—1 cc to 15 cc—of medication
- Less expensive and safer to use than processing and sterilizing a reusable system
- Fast and easy set up allows for immediate use
- Effortless delivery mechanism provides a controllable variation in the application rate of medication
- Directional tip provides healthcare workers with the ability to effortlessly direct and deliver medication to the desired location
- Its power sprayer feature blows the atomized medication down past the cords
- The atomized medication deeply penetrates the nasal cavity
- Non-stick trigger assembly overcomes problems common to spraying crystallizing, viscous medications.
- Quick and thorough anesthetization allows for instant scoping of patients



### Ideal Areas of Use

- Anesthesia
- Pulmonary Function/Respiratory Care
- GI Lab
- Endoscopy
- Radiation Oncology
- Critical Care
- Emergency Medicine



[AlcoveMedical.com](http://AlcoveMedical.com)

800.709.9888

HCPCS Code  
A9999  
Model EZ100m