

MoliCare[®]

Bed Mat & Rectangular

Cover all needs for
extra protection





MoliCare Bed Mat Eco

Disposable bed mat
with absorbent core made of
recycled cellulose layers

MoliCare Premium Bed Mat

Disposable bed mat
with absorbent core made of
cellulose fluff



Extra protection

absorbs liquid and provides extra
protection for various surfaces



Hygienic care

single-use secure protection



Dermatologically tested

dermatologically tested



Soft

with a soft non-woven top surface



Eco-friendly

absorbent layers are made of recycled cellulose



Pre-folded/compressed

for easy use and storage

Pre-folded



* 40 x 60 variants not compressed


The MoliCare Bed Mat Eco

Environmentally friendly, disposable bed mats
with an absorbent core of recycled cellulose layers,
providing extra leakage protection

Absorbs liquid and provides extra
protection for various surfaces



Absorbent layers
made of recycled
cellulose




Waterproof polyethylene
back sheet



MoliCare Bed Mat Eco 5 drops

The MoliCare Premium Bed Mat 5, 7 & 8 drops

Disposable bed mat with an absorbent core made of cellulose fluff providing hygienic care and extra protection.



MoliCare Premium Bed Mat 5, 7 & 8 drops

NEW

MoliCare Premium Bed Mat 9 drops

Newly developed, highly absorbent & soft disposable bed mat, providing hygienic care and comfort with extra leakage protection.

Absorbent core made of cellulose fluff and superabsorbent polymer

Superabsorbent core safely holds liquid and provides extra protection for various surfaces



**HYGIENIC
CARE**



Single-use
secure protection

Dermatologically
tested



Soft top sheet with
discreet, textile-like
back sheet



Textile-like back sheet prevents
slippage for secure hold in bed



MoliCare Premium Bed Mat 9 drops

The MoliCare Rectangular

MoliCare Rectangular is a hygienic solution for patients with faecal smearing or with light incontinence.


With their highly **absorbent core**, MoliCare Rectangular pads provide **secure protection** for everyone.

MoliCare Rectangular is **dermatologically tested**.

Soft
top sheet



Dermatologically
tested



Recommended
with
MoliCare
Fixpants



MoliCare Rectangular
3 drops

MoliCare Bed Mat & Rectangular advantages



* All Moli Care Rectangular pads and Moli Care Premium Bed Mat products

MoliCare® Bed Mat & Rectangular

Range overview

			USL Code	HARTMANN Code	Size	Packaging	Absorbency
MoliCare Premium Bed Mat		5 5 drops (Green Backing)	161065	161 065	60 x 90cm	120 (30 x 4)	850ml
		7 7 drops (Green Backing)	1206	161 066	40 x 60cm	150 (25 x 6)	561ml
			1207	161 070	60 x 90cm	100 (25 x 4)	1200ml
		8 8 drops (Green Backing)	8204	161 088	60 x 90cm	120 (30 x 4)	1395ml
	9 9 drops (White Backing)	8203A	161 074	40 x 60cm	90 (30 x 3)	1150ml	
MoliCare Bed Mat Eco		5 5 drops (Green Backing)	1205	159 001	40 x 60cm	300 (300 x 1)	490ml
MoliCare Rectangular		3 3 drops (Yellow Backing)	7165	160 890	20 x 40cm	300 (50 x 6)	440ml



For further information on MoliCare continence products, please contact your USL Account Manager or call customer services on **0800 658 814**



Going further
for health