



The Rocket Hysterocone™ Cannula has these new features:

- **MOVEABLE SILICONE CONE** ensures a soft, trauma free seal at the external OS.
- **1cm GRADUATIONS** on the distal 7cm to allow cone to be accurately positioned to match each patient's needs.
- **SMOOTH TIP & 6FG TUBE** ensures simple, low trauma insertion, ideal for the nulliparous cervix.
- **REMOVABLE HOLLOW STIFFENING TUBE** to aid insertion in and maintain rigidity during contrast injection.
- **CONTRAST PATH TAP** ensures retention of contrast when further volumes of media are required.
- **COMPATIBLE** with all aqueous contrast media.

Description	Size	Order Codes
Rocket® Hysterocone™ (HSG & Sonohystogram) Cannula Active length 280mm, with 13mm cone, removable hollow stiffener tube, contrast path tap. Sterile, for single use, peel pouch & tray packed in cartons of 10 cannula.	6FG	R57430-CONE USL: 857430-CONE

The Rocket® Hysterocone™ Cannula is a new concept in uterine contrast injection.

It combines the features of traditional cone and HSG catheters into a simple to use, comfortable and cost effective method of performing uterine canalisation for hysterosalpingography and hysterosalpingo-contrast sonography.

