



Package

- 20 punches in a dispenser box
- 50 boxes per shipping carton

Regular type

REF SIZE

BP-10F 1mm

BP-15F 1.5mm

BP-20F 2mm

BP-25F 2.5mm

BP-30F 3mm

BP-35F 3.5mm

BP-40F 4mm

BP-50F 5mm

BP-60F 6mm


BP-80F 8mm

Features

- Seamless stainless steel sharp cutting edge
- Ribbed handle for sure grip and control
- Size embossed on each punch for easy identification
- Wide range of sizes
- Sterile and ready to use
- Disposable

Clinical areas

- Dermatology
- Gynecology
- Podiatry



Long type
BP-L15K 1.5mm
BP-L20K 2mm
BP-L30K 3mm
BP-L35K 3.5mm
BP-L40K 4mm
BP-L50K 5mm
BP-L60K 6mm
BP-L80K 8mm





Features

- Stainless steel semi sharp cutting edge for soft or weak tissue
- Size embossed on each curette for easy identification
- Wide range of sizes
- Sterile and ready to use
- Disposable

Used in the following minor skin surgical procedures

- Curettage of a basal cell carcinoma
- Curettage of a verruca
- Curettage of other benign neoplasms of the skin

Clinical areas

- Dermatology
- Surgery
- Podiatry
- Veterinary

Package

- 20 curettes in a dispenser box
- 20 boxes per shipping carton

REF	SIZE
MK402	2 mm



MK403	3 mm
-------	------



MK404	4 mm
-------	------



MK405	5 mm
-------	------



MK407	7 mm
-------	------



A dermal curette biopsy technique for various dermatological diseases



A punch (Trepan) biopsy technique for various dermatological diseases

References : Manabu MAEDA, MD & Ph D *Department of Dermatology, Prefectural Gifu Hospital, 4-6-1 Noishiki, Gifu City, 500-8717, Japan*

Product specifications may be changed without notice



kai industries co.,ltd.
 1110,Oyana, Seki City, Gifu Pref. 501-3992 Japan
 Phone:+81 (0) 575 28 6600 Fax:+81 (0) 575 28 6611
 Email:info-ind@kai-group.com http://www.kai-ind.co.jp



kai Europe GmbH
 Löhndorfer Str. 51, D-42699 Solingen, Germany
 Phone:+49 (0) 212 232380 Fax:+49 (0) 212 2323899

DISPOSABLE BIOPSY PUNCH WITH PLUNGER SYSTEM




STERILE • STAINLESS STEEL • SINGLE USE • COLOR CODED • PATENT PENDING



Seamless Stainless Steel Cutting Edge
Ensures a consistently smooth cutting surface, allowing optimal specimen quality

Internal Plunger System

Ejects lodged skin specimen inside the metal lumen of the punch



	BPP-10F 1.0mm w/plunger	
	BPP-15F 1.5mm w/plunger	
	BPP-20F 2.0mm w/plunger	
	BPP-30F 3.0mm w/plunger	
	BPP-40F 4.0mm w/plunger	

Punches are available in three sizes, marked clearly on the handle, blister pack and box
Punches are individually packed in color-coded, transparent, sturdy blister packs, 20 per dispenser box

Product specifications may be changed without notice



CE
0197



kai industries co.,ltd.

1110, Oyana, Seki City, Gifu Pref. 501-3992 Japan
Phone : +81 (0) 575 28 6600 Fax : +81 (0) 575 28 6611
Email: info-ind@kai-group.com http://www.kai-ind.co.jp



kai Europe GmbH

Löhdorfer Str. 51, D-42699 Solingen, Germany
Phone : +49 (0) 212 232380 Fax : +49 (0) 212 2323899