

USL WRAPAROUND HINGED KNEE

The **USL Wraparound Hinged Knee** brace is designed to provide support, stabilisation and compression for most ACL, MCL, LCL strains, injuries and instabilities. Wraparound and hinged style accommodates all leg shape variations.

NEW



Open popliteal
eliminates bunching
and increases comfort

Features

- 1/8" neoprene construction.
- Open patella and open popliteal
- Front contact closures allow for hassle-free application
- Wraparound style accommodates all leg shape variations
- Hinges are removable if desired

Indications

- LCL and MCL strains, injuries and instabilities
- Meniscal injuries
- Mild OA



Find your size

USL Code	Size	Knee Circumference (15cm above knee centre)
23033	X-Small	38 - 44.5cm
23034	Small	44.5 - 49.5cm
23035	Medium	49.5 - 57.5cm
23036	Large	57.5 - 63.5cm
23037	X-Large	66 - 73.5cm
23038	XX-Large	73.5 - 81cm
23039	XXX-Large	81 - 89cm

