



U&L equipment

EQUIPMENT PRODUCT GUIDE



HUNTLEIGH CTG MONITORS



HUNTLEIGH CTG MONITORS

Sonicaid fetal monitors

The Huntleigh Sonicaid™ range of fetal monitors delivers innovative world-class obstetric solutions for midwives and obstetricians. Our extensive range allows monitoring throughout pregnancy in all environments.

Our success in fetal monitoring is born from more than 40 years' experience and whilst recognised as a major supplier of CTGs worldwide, we are continuously investing in new and improved solutions to meet the changing demands of our customers.

The range is complemented by our software systems which include maternity information, labour management and central review and archiving to satisfy all your hospital requirements.

**FOR MORE INFORMATION CONTACT:
USL MEDICAL CUSTOMER SERVICES ON 0800 658 814**



HUNTLEIGH SONICAID BD4000XS

Entry level monitoring

The BD4000xs series fetal monitors provide solutions for all areas of fetal monitoring. Using state of the art digital signal processing technology with enhanced FHR performance. The range includes antenatal and intrapartum models with a range of options and interface capabilities. This provides the most cost effective and flexible solution for all your fetal monitoring needs. Its compact size and styling belie its ruggedness and robust design, using the same materials, construction and design expertise as our world-renowned handheld dopplers.

Features:

- Twins capability as standard
- Unique no compromise full size twin FHR scale printing
- User adjustable, tachycardia, bradycardia, loss of contact and cross-channel rate verification alerts
- Hospital Information System and Vital Signs Monitor interface
- User-customisable trace annotation facility
- Automatic fetal movement detection
- Intuitive user interface and controls simplifying training
- Maternally sensed fetal movement marker
- FECG and IUP upgrade options



HUNTLEIGH SONICAID TEAM SERIES

Advanced fetal monitoring with CTG analysis

The Sonicaid Team range is designed for simplicity in fetal monitoring. Developed in collaboration with leading obstetric specialists, the small, easy to use and lightweight unit can be transported effortlessly - within hospitals, to clinics, doctors' offices or even to the home.

Features:

- Twin fetal heart rate scales provide the clearest differentiation between two heart rates and avoids the difficulty of interpreting offset or overlapping traces
- Obstetric management interface as standard
- Scrolling display of fetal heart rate and uterine activity traces
- Six hours of trace storage
- Records can be annotated at the touch of a button
- Maternally sensed and automatic (Actogram) fetal movement detection
- Tachycardia, bradycardia, signal loss and cross-channel rate verification alarms
- Twins (TeamCare Duo) and Intrapartum (Team IP) variants
- Care Antepartum CTG analysis provides an objective, numerical and consistent assessment of the CTG. Clinicians can rely on a large database of collective experience to support their own visual assessment
- Trend provides clear periodic changes over time in the FHR. For use during the first stage of labour, clinicians are quickly highlighted to changes in the variation, basal heart rate and deceleration size



HUNTLEIGH SONICAID FM800 ENCORE SERIES

High performance acute fetal/maternal intrapartum monitoring

The Sonicaid FM800 Encore is a high performance fetal and maternal monitor, using today's most advanced technology to produce the most accurate acquisition of data for clinician decision support. Whether you need standard fetal monitoring or equipment for specialist, high-risk care, the Sonicaid FM800 Encore series gives you accuracy and clarity in every situation.

In the fast-paced environment of the labour ward, you need equipment that's easy to use and data that is quick to read. The Sonicaid FM800 Encore gives you comprehensive fetal/maternal monitoring and analysis in one technologically advanced yet simple to use monitor.

- Using the very latest in ultrasound technology, the Sonicaid FM800 Encore gives increased performance. The wider, deeper, uniquely designed ultrasound beam shape makes it easy to locate and maintain contact with the fetal heart in mothers of all sizes
- Superbly accurate fetal ECG that works in the most demanding of environments - perfect for the busy labour room
- User defined alarm settings are held in memory for each new recording
- Unique Intelligent Fetal Analysis (IFA), providing fetal data analysis to support professional assessments
- Care Antepartum CTG analysis provides an objective, numerical and consistent assessment of the CTG. Clinicians can rely on a large database of collective experience to support their own visual assessment
- Trend intrapartum analysis provides clear periodic changes over time in the FHR. For use during the



PROUDLY NZ OWNED

USL has been built by creating seamless relationships. Since 1984 we've built enduring partnerships with customers, suppliers and our community. We set out to create the very best environment for our team so that they might deliver their unique skills to our markets.

WHAT YOU CAN EXPECT FROM US:

- We offer our customers dependability and experience they can rely on and trust.
- We deliver under IFOTIS: In-Full, On-Time, In-Specification.
- We operate a management reporting system that provides information on the frequency, recency, amount and type of purchases for all customers, and thereby enhancing their business purchasing efficiency.
- Normal guarantees apply and are rigorously honoured, and we guarantee to supply our key clients with critical stock items under IFOTIS.
- Our product range is comprehensive. It meets our Healthcare customers' needs, so they can confidently have a one supplier relationship for all medical products.
- Private Market
- General Practitioners
- Nurse Practitioners
- Private Hospitals (Medical, Surgical)
- Rest Homes
- Accident and Emergency Clinics
- Emergency Services
- Occupational Health Services
- Community Trusts
- Physiotherapists and Sports Medicine Clinics
- Sports Clubs and Organising Bodies
- Veterinary Services
- First Aid Suppliers
- Medical Laboratories
- Public Hospitals
- Specialist Clinics
- Community and Public Health Services
- Pharmacy Retail
- Industrial and Safety
- Dermatologists and Skin Clinics

OUR COMPANY AT A GLANCE:

- USL is a limited liability company, proudly New Zealand owned
- Our annual sales exceed \$NZ 60 Million
- Our head office and main distribution centre are based in Auckland
- The South Island is serviced from our Christchurch distribution centre
- We service 3,750 customers at all locations throughout New Zealand including all 20 of the National District Health Board



Order By Freephone:
Auckland: 0800 658 814
Christchurch: 0800 804 546

Order By Freefax:
Auckland: 0800 830 660
Christchurch: 0800 804 547

Auckland Head Office
494 Rosebank Road, Avondale, Auckland 1026
Email: customerservices@uslmedical.co.nz