




U&L equipment

EQUIPMENT PRODUCT GUIDE



ELECTROTHERAPY



INTELECT® MOBILE ULTRASOUND

- 1 and 3MHz frequencies
- Pulsed and continuous therapy operation (10%, 20%, 50% and 100%)
- Variable duty cycle outputs of 16Hz, 48Hz, or 100Hz
- Clinical protocols, the most appropriate application and parameter set for the patient's condition
- 10 user-defined memory positions
- Beam non-uniformity ratio (BNR) 6:1 maximum
- Lightweight design, along with battery-powered option
- Unit uniquely designed for tabletop, wall mount, or mobile use
- Ergonomically designed soundheads in variety of sizes (1cm², 2cm², 5cm² and 10cm²)
- All sound heads are interchangeable due to Electronic Signature™ residing on all transducers
- Head warming feature and head contact control with visual and auditory feedback on all sound heads
- 3 year warranty (accessories 1 year)

Standard Accessories

- 5cm² Applicator
- Ultrasound Gel
- User Manual on CD



INTELECT® MOBILE INTERFERENTIAL

- 2 channels of stimulation output
- Independent intensity and parameter controls for each channel
- Multiple waveforms on board: (2 Pole IFC, 4 Pole IFC, Galvanic, High Voltage Pulsed Current (HVPC), Microcurrent, Russian, Träbert, Monophasic Triangular, Monophasic Rectangular, Diadynamic, VMS and TENS)
- 15 user-defined memory positions
- Lightweight design, along with battery-powered option
- Unit uniquely designed for tabletop, wall mount, or mobile use
- 3 year warranty (accessories 1 year)

Standard Accessories

- Nylatex® Wrap.
- 6 x 8 cm carbon electrodes (x4)
- 6 x 8 cm sponges (x4)
- User Manual on CD



INTELECT® MOBILE ULTRASOUND

CODE	DESCRIPTION	UNIT
8558	Intellect® Mobile Ultrasound 1 and 3MHz	each

INTELECT® MOBILE INTERFERENTIAL

CODE	DESCRIPTION	UNIT
8560	Intellect® Mobile Interferential 2 Channels	each



INTELECT® MOBILE COMBINATION UNIT

- 1 and 3 MHz ultrasound frequencies
- Pulsed and continuous ultrasound therapy operation (10%, 20%, 50% and 100%)
- Variable duty cycle outputs of 16Hz, 48Hz, or 100Hz for ultrasound
- Multiple waveforms on board: (2 Pole IFC, 4 Pole IFC, Galvanic, High Voltage Pulsed Current (HVPC), Microcurrent, Russian, Träbert, Monophasic Triangular, Monophasic Rectangular, Diadynamic, VMS and TENS)
- 15 user-defined memory positions
- Lightweight design, along with battery powered option
- Unit uniquely designed for tabletop, wall mount or mobile use
- Ergonomically designed soundheads for ultrasound therapy in variety of sizes (1cm², 2cm², 5cm² and 10cm²)
- All sound heads are interchangeable due to Electronic Signature™ residing
- 3 year warranty (accessories 1 year)



*Therapy cart sold separately

INTELECT® MOBILE COMBINATION UNIT		
CODE	DESCRIPTION	UNIT
8561	Intellect® Mobile Combination Unit	each



INTELECT® ADVANCED COMBO THERAPYSYSTEM

Clinically convenient, the Intellect® Advanced Therapy System provides all the necessary tools in one package to maintain efficiency in the clinic and improve patient care.

- More than 25 clinical waveforms
- 2 independent Electrotherapy channels
- sEMG and sEMG + electrical stimulation standard
- High resolution colour display of all screens
- Over 200 clinical protocols
- Over 100 user-defined protocol slots
- Over 10 quick-link indications for quick parameter retrieval
- Impressive clinical library which includes full colour graphic anatomical and pathological libraries on vibrant LCD screen
- Documentation of treatment outcomes with patient data cards
- Sequencing to direct a designated therapeutic course without having to modify the unit
- 1 and 3Mhz ultrasound frequencies
- Pulsed and continuous modes (10%, 20%, 50% and 100%)
- Variable duty cycles (16Hz, 48Hz and 100 Hz)
- BNR 5:1 maximum
- Ergonomically designed sound heads in 1cm², 2cm², 5cm² and 10cm²
- 3 year warranty (accessories 1 year)

*Therapy cart sold separately

INTELECT® ADVANCED COMBO THERAPY SYSTEM		
CODE	DESCRIPTION	UNIT
9169	Intellect® Advanced Combo Therapy System	each



INTELECT® MOBILE BATTERY PACK

Electrotherapy on the go. A self-contained nickel metal hydride battery that stores conveniently in the bottom of the Intelect® Mobile Ultrasound, Stimulation, Combination, or Laser Unit.

INTELECT® MOBILE BATTERY PACK		
CODE	DESCRIPTION	UNIT
7113	Intelect® Mobile Battery Pack	each



INTELECT® ULTRASOUND APPLICATORS

Pre-calibrated, waterproof and ergonomically designed sound heads which operate at 1 and 3Mhz frequencies. The applicators are equipped with coupling indicators and are interchangeable between all Intelect® Mobile Units.

INTELECT® ULTRASOUND APPLICATORS		
CODE	DESCRIPTION	SIZE
9725	Intelect® Ultrasound Applicators	1cm
904781	Intelect® Ultrasound Applicators	2cm
9723	Intelect® Ultrasound Applicators	5cm



INTELECT® THERAPY CART SYSTEM

Designed with functionality in mind, the cart provides from 4 to 6 bins for storage (depending upon whether or not vacuum option is added). The cart integrates with the main unit to minimise cord and cable hassles. (Unit not included).

INTELECT® THERAPY CART SYSTEM		
CODE	DESCRIPTION	UNIT
8559	Intelect® Therapy Cart System	each



INTELECT® MOBILE CART ADAPTER

The mobile adapter allows the range of Intelect® Mobile Units to connect to the Intelect® Therapy Cart and Interferential suction unit.

*Indent Item

INTELECT® MOBILE CART ADAPTER		
CODE	DESCRIPTION	UNIT
2883*	Intelect® Mobile Cart Adapter	each



INTELECT® TENS STANDARD PORTABLE STIMULATOR

For symptomatic relief and management of chronic pain and acute post-surgical or post-trauma pain.

- Dual, independent channels
- Adjustable parameters
- 3 stimulation modes (burst, normal, modulation)
- Flip cover to protect top controls

INTELECT® TENS STANDARD PORTABLE STIMULATOR		
CODE	DESCRIPTION	UNIT
6948	Intellect® Tens Standard	each



INTELECT® NMES DIGITAL AND PORTABLE STIMULATOR

For the prevention of disuse atrophy as well as muscle re-education, increasing range of motion and increasing circulation.

- Neuromuscular electrical stimulation
- Dual, independent channels
- Adjustable parameters
- 3 stimulation modes (synchronous, constant, alternate)
- Patient compliance timer
- LCD display on digital model

INTELECT® NMES DIGITAL PORTABLE STIMULATOR		
CODE	DESCRIPTION	UNIT
9321	Intellect® NMES Digital Stimulator	each
10332	Intellect® NMES Portable Stimulator	each



HYDROCOLLATOR® E-1 STATIONARY HEATING UNIT

The benchmark against which all other heating units are judged.

- Thermostatically controlled temperature
- Dimensions (H x L x W) 41 x 33 x 20cm
- Includes 4 standard size HotPacs
- Main power - 110-240v 50/60Hz
- No plumbing required
- High quality stainless steel
- Weight - 7kg
- Tank capacity - 15 litres

HYDROCOLLATOR® E-1 STATIONARY HEATING UNIT		
CODE	DESCRIPTION	SIZE
3957	Hydrocollator® Heating Unit	15 litre



SONICLEAN HEAT PAN

The Soniclean Heat Pan is specifically designed as a water heater for the activation of thermoplastic materials. It comes complete with a temperature control and digital display.

SONICLEAN HEAT PAN	
CODE	DESCRIPTION
9584	Soniclean Heat Pan



VEINOPLUS SPORT®

Created in 2004 by Paris-based Ad Rem Technology, VEINOPLUS® device family has opened new horizons for the use of electro-stimulation in MEDICAL, SPORTS and WELL-BEING fields thanks to its innovative “Second Heart” technology. After 3 years of R & D and in co-operation with the top athletes, Ad Rem Technology launched for London Olympic Games a decisive innovation in sport recovery: VEINOPLUS® SPORT.

3 Effects

- Accelerated lactate clearance
- Local analgesic effect
- Relaxing action

2 Recoveries

- Whole body recovery: activates the Second Heart® (ie the calf-pump)
- Local recovery: relieves pain, prevents cramps and muscle soreness

1 Program

- Signal optimized for sport recovery
- Takes only a few minutes to be efficient

VEINOPLUS® SPORT

CODE	DESCRIPTION	SIZE
16003	Veinoplus Sport EMS Device	Each
16000	Veinoplus Device Electrodes	Each
16005	Veinoplus Sport Power Cuffs	Pair



VEINOPUS BACK®

The Veinoplus Back offers a unique program of pain relief and a relaxing massage which acts on the three main areas of the back. Compact and easy to use, you can take it anywhere. The signals of stimulation are a combination of TENS and NMES. TENS (Transcutaneous Electrical Nerve Stimulation) is widely recognised to alleviate pain whereas NMES (Neuro-muscular Electrical Stimulation) stimulates the muscles to cause their contractions and in turn relieve (relax) muscle spasms. Veinoplus® Back combines 4 effects providing back pain relief in one single stimulation program.

- Soothes back pain
- Delivers a relaxing massage
- Relieves muscle spasms
- Release natural endorphins to relieve general pain

VEINOPUS® BACK		
CODE	DESCRIPTION	SIZE
16001	Veinoplus Back EMS Device	Each
16000	Veinoplus Device Electrodes	Pair



PROUDLY NZ OWNED

USL has been built by creating seamless relationships. Since 1984 we've built enduring partnerships with customers, suppliers and our community. We set out to create the very best environment for our team so that they might deliver their unique skills to our markets.

WHAT YOU CAN EXPECT FROM US:

- We offer our customers dependability and experience they can rely on and trust.
- We deliver under IFOTIS: In-Full, On-Time, In-Specification.
- We operate a management reporting system that provides information on the frequency, recency, amount and type of purchases for all customers, and thereby enhancing their business purchasing efficiency.
- Normal guarantees apply and are rigorously honoured, and we guarantee to supply our key clients with critical stock items under IFOTIS.
- Our product range is comprehensive. It meets our Healthcare customers' needs, so they can confidently have a one supplier relationship for all medical products.
- Private Market
- General Practitioners
- Nurse Practitioners
- Private Hospitals (Medical, Surgical)
- Rest Homes
- Accident and Emergency Clinics
- Emergency Services
- Occupational Health Services
- Community Trusts
- Physiotherapists and Sports Medicine Clinics
- Sports Clubs and Organising Bodies
- Veterinary Services
- First Aid Suppliers
- Medical Laboratories
- Public Hospitals
- Specialist Clinics
- Community and Public Health Services
- Pharmacy Retail
- Industrial and Safety
- Dermatologists and Skin Clinics

OUR COMPANY AT A GLANCE:

- USL is a limited liability company, proudly New Zealand owned
- Our annual sales exceed \$NZ 60 Million
- Our head office and main distribution centre are based in Auckland
- The South Island is serviced from our Christchurch distribution centre
- We service 3,750 customers at all locations throughout New Zealand including all 20 of the National District Health Board



Order By Freephone:
Auckland: 0800 658 814
Christchurch: 0800 804 546

Order By Freefax:
Auckland: 0800 830 660
Christchurch: 0800 804 547

Auckland Head Office
494 Rosebank Road, Avondale, Auckland 1026
Email: customerservices@uslmedical.co.nz