



There is **no gentler way**
to close wounds

And now with extra strength
High Viscosity DERMABOND Mini
closing difficult wounds is even easier

DERMABOND Mini is now a High Viscosity formulation

- It is 6 times thicker than original **DERMABOND**
- Ideal for small wounds where sutures would otherwise be used
- It is easy to apply and, with 2 layers, provides fast closure
- It has all the benefits of **DERMABOND**
 - provides a microbial barrier¹
 - is water-resistant
 - offers flexibility and comfort for the patient



DERMABOND Adhesive Mini provides just the right amount of adhesive for most emergency room or general practice application

Code		Length of wound for 1 ampoule	Layers required
AHVM12	High Viscosity DERMABOND Mini	up to 4 cms	2

Fast precise closure, that also protects

High Viscosity **DERMABOND** Mini

- Eliminates or reduces the need for sutures
- Needs no refrigeration
- Provides the strength and security of healed tissue at 7 days
 - in less than 3 minutes²

**Please call customer service in Australia on 02 9878 9260
or 1800 252 194 or New Zealand 0800 803 988.**

1. Bhende S et al, In Vitro Assessment of Microbial Barrier Properties of DERMABOND Topical Skin Adhesive. Surgical Infections. 2002;3(3):251-257. 2. Date on File, Johnson & Johnson. * Trademark of Johnson & Johnson. Johnson & Johnson Medical Pty. Ltd. 1-5 Khartoum Road, North Ryde, NSW 2113. WH J113261

DERMABOND*
HIGH VISCOSITY
TOPICAL SKIN ADHESIVE 2-Octyl Cyanoacrylate
Every incision, every time.