

MATERIAL SAFETY DATA SHEET

SURGI gel™ PLUS STERILE LUBRICATING GEL

Page 1:3

1. Identification of the Product and Company

Product Name:	SURGI gel™ PLUS STERILE LUBRICATING GEL
Product Code:	SUR03038F, SUR02920F, SUR03029F
Other Names:	None allocated
Use:	To lubricate a body orifice to facilitate the entry of a surgical, diagnostic or therapeutic device.
Company Name & Contact Details	Medical Information Associate Distributed by Perrigo Australia Orion Laboratories Pty. Ltd., trading as Perrigo Australia ABN 56 009 293 136 25-29 Delawney Street, Balcatta, Western Australia 6021 AUSTRALIA Telephone (all hours): +618 9441 7800 FREE PHONE: 1800 805 546 FREE FAX: 1800 004 110 EMAIL: customerservice@perrigo.com.au ; WEBSITE: www.perrigo.com.au ORION [®] is a registered trademark of Orion Laboratories Pty Ltd
Other Information	All reasonable care has been taken to ensure information and advice contained in this data sheet is accurate at time of printing. However, Orion Laboratories Pty Ltd accepts no liability for any loss or damages suffered as a consequence of reliance on the information contained herein.

2. Hazards Identification

Hazard Classification	This product is not dangerous for the purposes of transport.
Risk phrase(s)	No component is present at sufficient concentration to require a hazardous classification This substance may cause irritation to eyes and may irritate skin of sensitive individuals.

3. Composition/Information on Ingredients

Chemical Entity	CAS No:	Proportion
Glycerol	56-81-5	10 - 30%
Carbomer 940	9003-01-4	< 10%
Propyl hydroxybenzoate	94-13-3	< 10%
Methyl hydroxybenzoate	99-76-3	< 10%
Purified Water	7732-18-5	60 - 100%

4. First Aid Measures

Inhalation:	Remove patient to fresh air. If breathing is difficult, apply artificial respiration. If required, seek medical advice.
Ingestion:	If swallowed and patient is conscious, wash out mouth with water. Seek medical advice and or contact a doctor or Poisons Information Centre (Australia 13 11 26).
Skin:	Wash off with water.
Eye:	If in eyes, hold eye lids open and flush with gently running water for at least 15 minute. If irritation occurs or persists, seek medical advice.
Advice to Doctor:	Treat symptomatically

5. Fire Fighting Measures

Extinguishing Media	Water spray, dry chemical powder, appropriate foam or carbon dioxide.
Hazards from Combustion products	Oxides of carbon (monoxide, dioxide) and nitrogen.
Precautions & Equipment for Fire Fighters	Wear self-contained breathing apparatus and protective clothing to prevent contact with skin and eyes.
Hazchem Code	None allocated.

MATERIAL SAFETY DATA SHEET

SURGI gel PLUS™ STERILE LUBRICATING GEL

Page 2:3

6. Accidental Release Measure

Personal Precautions	Wear appropriate clothing to avoid skin contact.
Clean Up of Spill	Mix with inert material (sand, vermiculite) and scoop into sealed containers for disposal. Area can be washed with soap and water.

7. Handling and Storage

Safe Handling Practices	Avoid contact with eyes. Keep away from ignition sources, heat, light and flame.
Storage	Store in well closed containers in a cool dry place out of direct sunlight.
Other Information	Incompatible with strong oxidizing agents and foods.

8. Exposure Controls; Personal Protection

Exposure Limits:	Not considered hazardous. There are no known Threshold Limit Values (TLV).
Engineering Controls	Precautions are not usually required
Personal Protection	Personal protection is not required.

9. Physical and Chemical Properties

Appearance & Odour:	Colourless transparent gel. Slight odour		
Flash Point:	> 100°C	Boiling Point:	Not known
pH:	6-7	Freezing/Melting Point:	Not known
Vapour Pressure:	Not known	Solubility:	Miscible with water
Vapour Density:	Not known	Specific Gravity or Density:	Not known

10. Chemical Stability and Reactivity Information

Conditions Contributing to Instability:	Stable under normal conditions of use.
Materials to Avoid:	Strong oxidizing agents.
Hazardous Polymerisation:	Will not occur.
Hazardous Decomposition Products:	Carbon monoxide, carbon dioxide, nitrogen oxides.

11. Toxicological Information

Inhalation:	No data available.
Ingestion:	No data available.
Skin:	May irritate skin of sensitive individuals.
Eye:	May cause irritation to eyes.

12. Ecological Information

MATERIAL SAFETY DATA SHEET
SURGI gel PLUS™ STERILE LUBRICATING GEL

Page 3:3

Mobility: No data available.

Persistence and Degradability: No data available.

Ecotoxicity: No data available.

13. Disposal Considerations

Disposal Methods & Containers: Empty containers should be decontaminated by rinsing with water prior to disposal.

Special Disposal for Landfill or Incineration: Waste material may be incinerated under controlled conditions where permitted. Refer to local Waste Management Authority Regulations for other approved methods.

14. Transport Information

UN Number: None allocated

UN Proper Shipping Name: None allocated

DG Class & Subsidiary Risk: None allocated

Packing Group: None allocated

Hazchem Code: None allocated

15. Regulatory Information

Registered as a device as per the Therapeutic Goods (Medical Devices) Regulations, 2002

16. Other Information

References:

- *Kept on file.*

END OF MSDS