



UNIVERSAL SHOULDER IMMOBILISER

- Polysling is a versatile shoulder immobiliser used to provide arm support or to immobilise the shoulder and arm



HOW POLYSLING WORKS

Polysling is used to provide arm support or to immobilise the shoulder and arm. Polysling is constructed from a strong, non-stretch material with a comfortable padded lining. The brushed nylon outer surface makes it possible to attach 'touch and close' fasteners as required, for maximum support and flexibility of use.

BENEFITS OF POLYSLING

- Versatile and effective sling material
- Reversible, for left or right sided use
- Only one size to stock – universal fitting
- Easily adjustable 'touch and close' fasteners
- Easy to apply and remove

MATERIAL PROPERTIES

- Polyurethane foam
- Brushed nylon cover
- Built-in 'touch and close' fasteners
- Natural rubber latex is not a constituent of Polysling or its packaging



AREAS OF USE

Polysling is a versatile and effective material which can be used in different configurations, such as left or right arm slings, or shoulder immobilisation slings.

POLYSLING ORDERING INFORMATION

Art. No.	Size	Pieces per inner	PIP Code	NHS Code
5555	Universal	6	TBC	TBC



Mölnlycke Health Care Ltd, Arenson Centre, Arenson Way, Dunstable, Bedfordshire LU5 5UL, UK. Tel. +44 [0] 870 6060766, www.molnlycke.co.uk
The Mölnlycke Health Care name and logo, Polysling™ are registered trademarks of Mölnlycke Health Care AB.
©2009 Mölnlycke Health Care. All rights reserved.

UK137550611