




The passive safety mechanism of VanishPoint® minimises the risk of needlestick injuries



BEFORE



AFTER

- Activation cannot be by-passed
- No change in clinical practice
- Needle automatically retracts into the syringe

When safety matters,
think TERUMO

VanishPoint®
Retractable Syringes

VanishPoint® Retractable Syringes

LATEX FREE

U-100 Insulin Syringes

Product Code	Colour	Volume	Gauge	Length		Carton
		mL	G	inches	mm	Qty
RTI*15291		0.5	30	$\frac{2}{10}$	5	800
RTI*15271		0.5	30	$\frac{5}{16}$	8	800
RTI*15221		0.5	30	$\frac{1}{2}$	12.7	800
RTI*10291		1	30	$\frac{2}{10}$	5	800
RTI*10271		1	30	$\frac{5}{16}$	8	800
RTI*10211		1	29	$\frac{1}{2}$	12.7	800

Shelf Pack = 100

Tuberculin Syringes

Product Code	Colour	Volume	Gauge	Length		Carton
		mL	G	inches	mm	Qty
RTI*10151		1	25	$\frac{5}{8}$	16	800
RTI*10161		1	25	1	25	800
RTI*10131		1	27	$\frac{1}{2}$	12.7	800

Shelf Pack = 100

Hypodermic Syringes

Product Code	Colour	Volume	Gauge	Length		Carton
		mL	G	inches	mm	Qty
RTI*10371		3	20	1	25	600
RTI*10381		3	20	$1\frac{1}{2}$	38	600
RTI*10361		3	21	$1\frac{1}{2}$	38	600
RTI*10311		3	23	1	25	600
RTI*10321		3	23	$1\frac{1}{2}$	38	600
RTI*10301		3	25	$\frac{5}{8}$	16	600
RTI*10391		3	25	1	25	600
RTI*13011		3	25	$1\frac{1}{2}$	38	600
RTI*10571		5	20	1	25	600
RTI*10581		5	20	$1\frac{1}{2}$	38	600
RTI*10551		5	21	1	25	600
RTI*10531		5	22	1	25	600
RTI*10541		5	22	$1\frac{1}{2}$	38	600
RTI*11041		10	22	$1\frac{1}{2}$	38	600

Shelf Pack = 100

Terumo Australia Pty Ltd

Macquarie Park NSW 2113 Australia
 Mount Albert Auckland 1025 New Zealand
 E: cs_australia@terumo.co.jp
 Orders: au_orders@terumo.co.jp

T: 1800 837 866
 T: 0800 66 77 57

VanishPoint® is a Registered Trademark of Retractable Technologies Inc.
 © 2017 Terumo Australia Pty Ltd. GHP170413D659V3 | PBD00208



Engineered for protection