



The PLUS Seat is designed to relieve muscle tension and reduce the body's natural pressure points for total sitting comfort.

By improving posture, increasing blood flow, and reducing pressure, you can be healthier when you sit. It works on any type of seat, from car seats to hard seats. Its contoured design relieves pressure on the hips so they remain upright in proper alignment thus reducing stress on the lower back. The PLUS Seat, made from waterproof closed-cell EVA material won't absorb moisture or generate heat for easy cleaning. Its patented design reduces seat surface temperatures and relative humidity by providing continuous airflow for ultimate comfort.

**AIRFLOW**

Port and pillar design provides continuous airflow to reduce temperature and moisture to the skin.

**WEIGHT DISTRIBUTION**

Bottom surface uses shorter pillars allowing the pelvis to immerse into the seat, alleviating peak pressure points.

**CORRECT POSTURE**

Contoured design helps keep the pelvis in proper alignment, reducing stress on the lower back.

**WATERPROOF**

Durable closed cell EVA is completely waterproof and heat resistant to 158F / 70C.