


WelchAllyn®



Welch Allyn
SureTemp® Plus 692 Electronic Thermometer



The SureTemp Plus 692 Thermometer— with the speed and accuracy you require, and the safety and security you demand.

1 Easy:

Intuitive user interface, easy-to-read LCD screen and flashing icons, helps select and identify temperature measurements and modes.

2 Safe:

Designed with interchangeable oral/axillary (blue) and rectal (red) probe wells to address patient safety and help reduce the risk of cross-contamination.

3 Fast and Accurate:

With proven accuracy, the SureTemp 692 thermometer captures quick readings—approximately 4-6 second oral, 10-12 second paediatric axillary and rectal and 13-15 second adult axillary measurements.

4 Secure:

To top it all off, our security features make sure the unit stays where it is supposed to be. A location ID on the display and other innovative features help prevent theft.

- Easy-to-use, intuitive operation
- Interchangeable oral/axillary and rectal probe wells help reduce the risk of cross-contamination



SureTemp Plus 692—setting the standard of care

ID field for electronic labeling of thermometer



User-selectable icon and/or words for identification of thermometer mode

Battery life indicator

Interchangeable probe with push-button probe cover ejection

Removable probe well for easy cleaning

C/F button for toggling between Fahrenheit and Celsius

Prominent recall button for display of last temperature taken

Ergonomic design with rubberised grips

Convenient storage for 25 probe covers

Large, easy-to-read, backlit LCD

Mode button for selecting oral, paediatric axillary or adult axillary

Pulse timer button/monitor mode

NEW feature makes it easier to switch to Monitor Mode for measuring a patient's thermal steady state (approx. 3 to 5 minutes)

Interchangeable rectal probe and well reduce risk of cross-contamination



The **SureTemp Plus 692** with **standard wall mount** has storage for an additional 50 probe covers.



The SureTemp Plus 692—setting the standard of care even higher

SureTemp Plus 692 Product Specifications

Calibration Accuracy	±0.1 °C Monitor Mode	
Temperature Predict Time	Oral: Adult Axillary (18 years and older): Paediatric Axillary (17 years and younger): Rectal:	Approximately 4-6 secs Approximately 15 secs Approximately 10 secs Approximately 10 secs
Unit Dimensions	215 mm x 81 mm x 62 mm	
Weight	357 g	
Ambient Temperature Range	10 °C to 40.0 °C	
Patient Temperature Range	26.7 °C to 43.3 °C	
Power Source	(3) 1.5Vdc AA batteries	
Battery Life	Approximately 5,000 temperatures	

Ordering Information:

THERMOMETERS

- 01692-400** SureTemp Plus 692 Electronic Thermometer, Wall Mount, Security System with ID Location Field, Oral Probe and Oral Probe Well with 1.2 m cord
- 01692-201** SureTemp Plus 692 Electronic Thermometer, Wall Mount, Security System with ID Location Field, Rectal Probe and Rectal Probe Well with 1.2 m cord
- 01692-500** SureTemp Plus 692 Electronic Thermometer, Wall Mount, Security System with ID Location Field, Oral Probe and Oral Probe Well with 2.7 m cord

ACCESSORIES

- 21326-0000** Wall Holder
- 21330-0000** Holder for Spare Probe and Well
- 06138-000** Calibration Key Assembly, 690/692
- 02892-000** Probe and Well Kit, 1.2 m Rectal
- 02892-100** Probe and Well Kit, 2.7 m Rectal
- 02893-000** Probe and Well Kit, 1.2 m Oral
- 02893-100** Probe and Well Kit, 2.7 m Oral
- 02891-0000** Probe Well Blue
- 02891-1000** Probe Well Red
- 08273-000** Thermometer Stand with Holder for Spare Probe and Well
- 21333-0000** Bracket, M600 Stand Upgrade, with Holder for Spare Probe and Well

REPLACEMENT PROBE COVERS

- 05031-101** Case of 1,000 Probe Covers
- 05031-105** Case of 5,000 Probe Covers

To find out more, please contact your **Welch Allyn Sales Representative** or Customer Service on **1800 650 083**



Welch Allyn Australia Pty. Ltd.
 Suite 4.01, 2-4 Lyonpark Road
 Macquarie Park NSW 2113 AUSTRALIA