



DONJOY[®]

Performance & Lifestyle

Stabilizing Pro Stabilizing Speed Pro

Ankle Braces

Designed to meet the needs of demanding athletes

The Stabilizing Pro and new Stabilizing Speed Pro ankle braces can reduce the risks of ankle injury by 61%¹. Both feature a new contoured design made of durable ballistic nylon, reinforced eyelets, a wide, circumferential compression strap and removable stays.

The Stabilizing Pro braces gives athletes the confidence that their ankles are protected.

Stabilizing Speed Pro



- Also available with new "speed-lacing" system that makes applying the brace exceptionally fast and easy
- Stabilizing Speed Pro features all the same benefits as the Stabilizing Pro version

¹ McGuire T.A., Brooks A. AM J Sports Med.2011 "The Effect of Lace-up Ankle Braces on Injury Rates in High School Football Players"

Stabilizing Pro

Features & Benefits

New Improved Low-Profile Design

- Fits easily into all shoe types including sports shoes

Wide (75 mm) Circumferential Elastic Strap

- Provides tibia/fibula compression and stability
- Secures all strapping

Strapping that performs consistently

- Non-stretch figure-8 straps lock calcaneus providing ideal inversion/eversion control

Padded, Mesh Tongue

- Provides comfort and helps to conform brace to foot

Indications

- Chronic ankle instabilities and/or sprains
- Prophylactic use

Universal Fit

- Fits left and right foot

Bilateral Stays (Inside)

- Removable medial/lateral stays provide additional support and stability without compromising comfort

Reinforced Eyelet area

- Helps prevent lace tearing

Non-stretch Ballistic 800D Nylon

- Robust material for tough athletes
- Delivers comfort, performance and durability






Stabilizing Pro / Stabilizing Speed Pro

ORDERING INFORMATION

STABILIZING PRO ANKLE		
PART NUMBER	COLOUR	SIZE
11-3234-x-06000	Black	XXS-XXXL
11-3234-x-15000	White	XXS-XXXL

STABILIZING SPEED PRO ANKLE		
PART NUMBER	COLOUR	SIZE
11-3235-x-06000	Black	XXS-XXXL

SIZES AVAILABLE (Measurement is circumference, taken from the base of the heel around ankle, while wearing an athletic sock)		
x = 0 = XXS	23-25 cm (9" - 10")	
x = 1 = XS	25-28 cm (10" - 11")	
x = 2 = S	28-30 cm (11" - 12")	
x = 3 = M	30-33 cm (12" - 13")	
x = 4 = L	33-36 cm (13" - 14")	
x = 5 = XL	36-38 cm (14" - 15")	
x = 6 = XXL	38-41 cm (15" - 16")	
x = 7 = XXXL	41-43 cm (16" - 17")	

DJO Global | 1a Guildford Business Park | Guildford | Surrey | GU2 8XG | UNITED KINGDOM
www.DJOglobal.eu



Together in Motion™