## TAMPER EVIDENT

# Rainbow Trays



**SECURE • AUDITABLE • ACCESSIBLE** 

# Rainbow Trays

The Tamper Evident Rainbow Tray<sup>™</sup> provides secure and auditable storage for all controlled drugs whilst in transit between secure locations.

As part of the Rainbow Tray<sup>™</sup> system, the Tamper Evident tray delivers the same organisational benefits as the Rainbow Trays<sup>™</sup> in addition to providing a robust tamper evident sealed unit that encapsulates the Chain-of-Custody information required by Accountable Offcers for Controlled Substances.



### TAMPER EVIDENT RAINBOW TRAYS™

### **SECURE**

- Securely sealed immediately after dispensing
- No access to contents after securing with label until required
- Crush resistant tray protects contents
- Comprehensive tamper evident features show any attempt to access the tray and its contents

### **AUDITABLE**

- Chain of Custody compliant
- Ideal for securing against diversion
- Label notes are sealed and cannot be altered
- Origination, destination and contents are clear and unambiguous
- Versatile write what you need on the label, make it your own.
- Date and time of dispensing and expiry with second check option

### **ACCESSIBLE**

- Contents are easy to identify in the transparent tray
- Easy removal from the base using tab
- Quick release perforation to open tray and access contents

### TAMPER EVIDENT RAINBOW TRAYS™ ARE COMPLIANT WITH

### **NHS England**

Controlled Drugs Regulations 2013

### **CQC - Care Quality Commission**

Guidelines on the Safer Management of Controlled Drugs 2016

# NICE - National Institue for Health and Care Excellence

NICE guideline [NG46] Published date: April 2016

### **RCoA - Royal College of Anaesthetists**

Guidelines on the Storage of Drugs in Anaesthetic Rooms 2016

# EBA - European Board of Anaesthesiology (EBA UEMS)

Recommendations for Safe Medication Practice



# TAMPER EVIDENT RAINBOW TRAYS™ SUPPORT



#### **ACCOUNTABLE OFFICERS**

The simple auditable label can be customised for local protocols to verify the chain of custody for the controlled drugs whilst in transit between Pharmacy and Patient. The tray and label are compliant with CQC & NHS recommendations for the security of controlled drugs.



#### **PHARMACISTS**

Easy to prepare and store, the trays and labels will verify the contents, origination and destination while referencing the original prescription. Ideal for dispensing for multiple locations and using the tamper evident label to complete the audit trail as required.



#### ANAESTHETISTS AND NURSES

The secure tray design allows the anaesthetist to confirm the authenticity of the tray contents and access them in a timely manner when required. Suitable for use where expired or unused drugs need to be safely returned to the pharmacy.



#### **EMERGENCY DEPARTMENT**

Ideal for use in Emergency Departments where drugs can be easily identified and accessed whilst remaining secure. The cutomisable label makes the tray suitable for drugs used for trauma. Customisable label.
Complete with
information as required.
Add text or stickers and
encapsulate for total
security.

Can also be used to return unused or expired drugs safely to pharmacy.

Ideal for Emergency Drugs where secure storage is required. Tamper Evident
Rainbow Trays™ are the
natural extention of the
Rainbow Trays™ system
to provide secure and safe
storage of Anaesthetic drugs
prior to administration whilst
supporting the critical care
colour labelling
system.

For infection control.
Bacteriostatic material will not harbour bacterial growth.

The Versatile design of the tray and the label allows the trays to be used for any drug type required and in any department.

Allows the pharmacy to draw up anaesthetic drugs in a controlled environment and distribute to operating theatres prior to procedures.

Quick and easy to open to access the drugs in an emergency.

The Tamper
Evident system allows
anaesthetists and nurses
to draw up their drugs in
advance and secure them
in the Tamper Evident
Rainbow Trays™ for
future use.

Tamper Evident security features:

- Permanent adhesive
  - Security cuts to evidence tampering
  - Encapsulated text

### WHAT WILL YOU USE YOURS FOR?

**CONTROLLED DRUGS** 

**OPIOIDS** 

**AMPOULES** 

**ANAESTHESIA** 

**OBSTETRICS** 

**CRITICAL CARE** 

**TRAUMA** 

PAIN MANAGEMENT

**SYRINGES** 

**GENERAL MEDICINES** 

**PHARMACY** 

**DAY CASES** 

LAB SAMPLES











