



**Rocket® HSG Catheters provide safe, effective and well-tolerated uterine canalisation for hysterosalpingography and hysterosalpingo-contrast sonography.**

#### CONTENTS:

- Rocket® HSG Catheter
- Introducer sheath
- Pre-attached 2ml syringe
- Protective balloon cover

#### Rocket® HSG Catheters have these well-established and trusted features:

- **LATEX FREE** balloon minimises risk of patient tissue reactions.
- **SMOOTH TIP & LOW PROFILE BALLOON** ensures low trauma insertion and removal – ideal for the nulliparous cervix.
- **REMOVABLE STIFFENING SHEATH** to aid insertion in the heavily retroverted or anteverted cervix.
- **CONTRAST PATH TAP** ensures retention of contrast when further volumes of media are required.
- **INTEGRAL WIRE STYLET** for improved malleability without using tiresome removable stylets (ST versions)
  - ◇ Improves ease of insertion
  - ◇ Increased value for money
- **COMPATIBLE** with all aqueous based contrast media.
- **COMPLETE KIT** for ease, speed and cost effective use of resources.

Description	Size	Order Code
<b>Rocket® HSG Catheter</b> HSG Catheter Length 280mm, Maximum balloon volume: 3ml. Sterile, for single use, in cartons of 10 units.	5FG	<b>R57422</b>
	7FG	<b>R57424</b>
<b>Rocket® HSG Catheter with integral wire stylet.</b> HSG Catheter length 280mm, Maximum balloon volume: 3ml. Sterile, for single use, in cartons of 10 units	5FG	<b>R57422-00-ST</b>
	7FG	<b>R57424-00-ST</b>