



The **Rocket® Seldinger Chest Drainage Set** is indicated for the percutaneous insertion of a pleural drain for fluid and air and is particularly suitable for: ¹

- Pneumothorax, spontaneous or iatrogenic
- Effusions, benign or malignant, except those with high protein content.
- Low viscosity empyema

The **Rocket® Seldinger Chest Drainage Set** features:

- Improved patient comfort using small bore catheterisation ^{2,3}
- Reliable insertion via Seldinger approach under direct vision, x-ray screening or ultrasound guidance
- 10cm dilators with safety stops to minimise the risks of over-insertion associated with the traditional 20cm dilator.
- Drapes and a gown are included in all PK coded packs in response to NPSA/2008/RRR3

References:

1. *Clinical Evaluation of Seldinger Drainage Kit.* Beaumont A C et al Proceedings of Intensive Care Society Meeting. December 1999.
2. *Catheters smaller than 24 French gauge can be used for chest drains.* Taylor PM. BMJ (1997) 315:186
3. *Management of Malignant Pleural Effusions* Antunes G, Neville E, Thorax (2000) 55:981-983

Description	Size	Code
Rocket® Seldinger Chest Drainage Set with Safety Guards Safety scalpel, echogenic tipped insertion needle, guidewire, 10cm dilators with Safety Guards, 30cm catheter with radiopaque stripe and tip, syringe, stopcock, fir-tree & wide-bore luer connectors. Sterile, single use in cartons of 5 units.	12FG	R54544-12-SG USL: 854544-12-SG
	18FG	R54544-18-SG USL: 854544-18-SG
Rocket® Seldinger Chest Drainage Tray Complete Ward Procedure Set Safety scalpel, echogenic tipped insertion needle, guidewire, 10cm dilators with Safety Guards, 30cm catheter with radiopaque stripe and tip, syringe, stopcock, fir-tree & wide-bore luer connectors, skin fixation set, ward procedure pack including drapes, swabs, sponges, 10ml syringes, gallipot. Sterile, for single use in cartons of 5 units.	12FG	R54544-12-PK USL: 854544-12-PK
	16FG	R54549-16-PK USL: 854549-16-PK
	18FG	R54544-18-PK USL: 854544-18-PK
	20FG	R54549-20-PK USL: 854549-20-PK
Rocket® Seldinger Chest Drainage Training DVD A clinician's guide to the safe insertion and connection of Rocket® Seldinger Chest Drainage Sets		R54544-DVD