



Welch Allyn CP 150™ Electrocardiograph



Simple. Fast. Connected.

The user-friendly CP 150 ECG—designed to help improve your facility workflow.

- > 7" Color, touch-screen display for fast and easy entry of accurate patient information
- > Instant-on feature powers the device up quickly and lets you take an ECG at a moment's notice
- > Advanced filters for optimal ECG trace quality ensure accurate readings
- > Three (3) lead ECG preview
- > Optional adult and pediatric MEANS ECG interpretation software offers a second opinion
- > Battery operation lets you easily transport your device to the patient
- > Store up to 100 test results on the device or transfer to a USB memory stick



TOUCH-SCREEN KEYBOARD

No need to track between the display and keyboard, saving time and helping reduce errors.



LEAD QUALITY DISPLAY

Allows you to quickly and easily identify lead connection problems.

WelchAllyn®




Welch Allyn CardioPerfect Workstation

Bring your cardiopulmonary data together with the Welch Allyn CardioPerfect® Workstation

- Manage patient information and tests performed with the full suite of Welch Allyn cardiopulmonary devices
- Get paperless efficiency and connectivity to your hospital information system (HIS) via HL7 interface
- Store, access, review, compare and edit diagnostic test data from any computer on your network
- Retrieve records with the touch of a key, and review and edit right on the screen



Flexible Connectivity Options



Connect, manage and share patient information

With a USB or Ethernet connection, you can transfer patient data from the CP 150 to the Welch Allyn CardioPerfect Workstation, giving you the ability to connect, save and manage patient data in your HIS.



Telemedicine

Enhance your diagnostic capabilities with telemedicine

With the CP 150, your hospital can provide support to remote recording sites requiring an expert opinion. Simply use the CardioPerfect Workstation to receive diagnostic testing data from multiple locations and get the capability to provide expert opinions and centrally manage patient information.



Multiple USB ports for connectivity to:

- > CardioPerfect software
- > Barcode reader
- > External memory
- > External printer

CP 150 ECG Mobile Cart

- > Attractive, ergonomic design
- > Hospital-grade wheels also available
- > Optional patient cable arm and shelf



Model Shown here with hospital-grade wheel option

Accessory Choices: AHA, IEC, reusable/disposable electrodes

Language: English, French, Italian, German, Spanish, Portuguese, Polish, Danish, Dutch, Chinese, Russian, Norwegian, Swedish, Vietnamese

Power Source: Available with plug styles that are compatible with any power source.

Contact your local Welch Allyn representative or visit www.welchallyn.com/CP150 today for more information.

Welch Allyn Corporate Headquarters
 4341 State Street Road, P.O. Box 220
 Skaneateles Falls, NY 13153-0220 U.S.A.
 (p) +1.800.535.6663 | (f) +1.315.685.0054

Welch Allyn European Customer Service: +353 46 90 67790

Australia: +61 29 638 3000
 China: +86 21 6327 9631
 France: +33 155 69 58 49
 Germany: +49 6950 985 132
 Latin America (Miami office): +1 305 669 9003

Malaysia: +603 78753341
 Netherlands: +31 202 061 360
 Singapore: +65 6419 8100
 Sweden: +46 858 536 551
 United Kingdom: +44 207 365 6780

WelchAllyn®
 Advancing Frontline Care™