

MATERIAL SAFETY DATA SHEET

Clinidet Medical & Dental Instrument & Equipment Detergent

Quality Medical Innovations Pty Ltd.

Clinidet is manufactured for Quality Medical Innovations Pty Ltd
Distributed By: Majac Medical Products Pty Ltd
PO Box 74 Virginia,
Queensland, Australia 4014
Tel: 07 3265 6355 (+61 7 3265 6355) Fax: 07 3865 2729 (+61 7 3865 2729)

1.0 IDENTIFICATION

PRODUCT NAME: **Clinidet** MEDICAL and DENTAL INSTRUMENT & EQUIPMENT DETERGENT

OTHER NAMES: None

USE: **Clinidet** is effective for the cleaning of instruments by manual or machine operations. General hard surface cleaning (benches and tables etc.), laboratory equipment and glassware.

2.0 HAZARDS IDENTIFICATION

DANGEROUS GOODS RATING: Not classified as hazardous according to criteria of Worksafe Australia

HAZCHEM CODE: None allocated **UN NO.:** None allocated

Dangerous Goods Class: None allocated

Poisons schedule: Not scheduled

3.0 INGREDIENTS

Ingredients not classified as hazardous according to criteria of Worksafe Australia

Non Ionic Surfactants	68551-13-3	<20%
Sodium Citrate	68-04-2	<5%
Sodium Metasilicate	10213-79-3	<5%
Triethanolamine	102-71-6	<5%

4.0 FIRST AID

Swallowed: Wash mouth with water and then give plenty of fluids to drink. Seek medical advice or phone Poisons Information on 13 11 26.

Eye: Hold eyes open and wash with water for at least 10 minutes. Remove contact lenses if worn. In all cases of eye contamination it is a sensible precaution to seek medical advice.

Skin: Rinse well with cold water.

Inhaled: Move patient to fresh air and seek medical advice if necessary.

Advice to Doctors: Treat symptomatically. For further information phone Poisons Information on 13 11 26.

5.0 FIRE FIGHTING MEASURES

Flammability: Not Flammable

MATERIAL SAFETY DATA SHEET

Clinidet Medical & Dental Instrument & Equipment Detergent

13.0 DISPOSAL CONSIDERATIONS

Detergent concentrate may be disposed by diluting with water then through the sewage system.

14.0 TRANSPORT INFORMATION

No specific requirements

15.0 REGULATORY INFORMATION

None allocated

16.0 OTHER INFORMATION

DISCLAIMER

This information is based on data believed to be accurate at the time of writing, but subject to change without notice. No warranty is expressed or implied as to the information's accuracy or completeness.