



CAREWELL  
**AGED CARE BED  
AND MATTRESS GUIDE**

ODEN AGED CARE BED

Are you designing or operating a home that has residents who are high care? Are they in need of a low low bed for minimal restraints, but height to allow your carers to provide care without risk of damaging their backs? Knee breaks and adjustable foot panels to manage blood pooling? Is wall damage costing your home more than another staff member? And do you have a limited budget in which to do all this?

If any of these questions apply to you then the feature packed Oden may be the answer for your home.

- X lift functionality
- Head tilt functionality
- Electric three function bed
- Electric knee break
- Central Locking 125 mm Castors
- Epoxy powder coating
- Safe Working Load 215 kg
- In built leg extension
- Wood grain Head and foot board
- Height range: 210 to 760mm
- Removable metal infection pads




ODEN AGED CARE BED	
CODE	DESCRIPTION
13651	3 section bed with central locking
NIMIR CARE BED	
14055	King single bed, 3 section
CROFT CARE BED	
13692	2 section bed with individual locking castors

Accessories:

14007	Carewell Side Rails
14224	Carewell Overhead Table
TBA	Carewell Self Help Pole


## ODEN TECHNICAL DETAILS



### Technical Information

Max Dimensions		Length	1960 - 2260 mm
		Width	910 mm
Patient surface dimensions			190-215 x 91cm
Backrest	- length		76 cm
	- tilting angle		85°
Headrest	- tilting angle		35°
Pelvic	- length		29 cm
Thigh rest	- length		30 cm
Calf rest	- length		65 cm
Height Adjustment	Lowest Height (bottom mattress platform)		21 cm
	Lowest Height (top mattress platform)		26 cm
	Top Height (top mattress platform)		76 cm
Safe Working Load			215 kg
Measured noise level			max 47dB
Electrical Details			
Supply Voltage			24V DC
Max. Load	- push		6000 N
	- pull		4000 N
Speed			4.2 mm/Sec
Duty Cycle			10% Max 2min/18 min

## Nimir Technical Details



## Technical Information

Max Dimensions	Length	1960 ~ 2260 mm
	Width	1060 mm
Patient surface dimensions		
Backrest	- length	190~215 x 91cm
	- tilting angle	85°
Headrest	- tilting angle	35°
Pelvic	- length	29 cm
Thigh rest	- length	30 cm
Calf rest	- length	65 cm
Height Adjustment	Lowest Height (bottom mattress platform)	35 cm
	Lowest Height (top mattress platform)	42 cm
	Top Height (top mattress platform)	76 cm
Safe Working Load		215 kg or 350kg
Measured noise level		max 47dB

## Electrical Details

Supply Voltage		24V DC
Max. Load	- push	6000 N
	- pull	4000 N
Speed		4.2 mm/Sec
Duty Cycle		10% Max 2min/18 min



### Pro-Flex Memory Foam Pressure Relieving Mattress

Designed to exceed even the most rigorous requirements of the acute care environment.

- Minimized shearing
- Outstanding comfort for longer , more restful sleep
- Superior in-bed mobility
- Waterlow Score 16.8
- Non turn mattress
- Non slip base
- Heel support
- Multi layered foam construction
- Mattress profiling
- Waterproof, breathable cover
- Stain resistant Fire resistant Antimicrobial
- Dust mite resistance
- Concealed Zip and Welded seams

### ProVex Concave Memory Foam

The CareWell convex mattress is a 195cm long and 90cm wide mattress with wings covered in a waterproof PU cover. Ideal for use with patients prefer not to have side rails.

- Polyurethane two way stretch fusion cover
- Non slip PVC bottom
- Waterfall zipper 1800 to prevent moisture
- Slow recovery visco-elastic memory foam for body contour
- Bottom layer: high density foam for body support
- Raised wing edging
- Heel support
- Mattress profiling
- No turning required
- Depth: 140 mm Length: 1950 Width 90 Wing 190 mm
- Fire rated

PRO-FLEX	
CODE	DESCRIPTION
13650	Pro Flex Memory Foam Mattress - Single
14054	Pro Flex Memory Foam Mattress - King single

PROVEX CONCAVE	
CODE	DESCRIPTION
13652	Pro Vex Concave Memory Foam Mattress - Single

# The CareWell Guarantee

Our goods come with a guarantee that you are entitled to a replacement or refund for a major failure and for compensation for any other reasonably foreseeable loss or damage. You are also entitled to have the goods repaired or replaced if the goods fail to be of acceptable quality and the failure does not amount to a major fault.

Key points:

- Three year warranty on the frame
- Three year warranty on all other parts including electronics
- Hand Control Excluded

CareWell manufactured parts carry a three year warranty against faulty workmanship, commencing from the date of manufacture. Parts not manufactured by CareWell carry the respective manufacturer's warranty with a minimum of three years, excluding batteries. The warranty includes the replacement of any parts that become defective during the warranty period through normal recommended use.

The warranty excludes repairs required arising out of any of the following – normal wear and tear, and abuse of product.







John Paul  
Product Manager

m: +64 21 728 359  
ddi: +64 9 829 0595  
e: [john.paul@uslmedical.co.nz](mailto:john.paul@uslmedical.co.nz)