




OSCMED
+ SWISS MEDICAL COMPANY



CONTENT

Cast cutter saw OSCIMED 2000.....	3
Cast cutter saw OSCIMED ERGO II.....	4
Vacuum cleaner for OSCIMED ERGO II and OSCIMED 2000	5
Cast cutter saw OSCIMED PSV	6
Vacuum cleaner for OSCIMED PSV	7
Plaster table and options	8
Plaster trolley.....	9
Plaster instruments.....	10
Knee rest for casting	
Sterilization baskets.....	11

OSCIMED 2000

Cast cutter saw




ref. OSC 20101

LIGHT WEIGHT AND LOW NOISE, PERFECT FOR CHILDREN PLASTER ROOM OR HAND SURGEONS

The OSCIMED 2000 cast cutter saws is equipped with a brushless motor for a precise and low noise use.

Rated power	100 Watts
Length of the saw body	210 mm
Effective diameter	46 mm
Weight without power supply	760 gr.
Noise level	62 dBA
Length of power supply cable	1,00 m
Length of the cable 24 VDC	3,50 m




Adjustable guide and depth gauge.



Usable with one hand thanks to its very short (210 mm) and lightweight (860 g) shape.



Very low noise level thanks to a patented anti-noise technology.



Blade change without tools.

SAW BLADES AND ACCESSORIES



Blade lock nut
ref. OSC 202



Vacuum nozzle
ref. OSC 240

Saw Blades	Ø 77 mm	Ø 65 mm	Ø 65 mm	Ø 50 mm	Ø 50 mm
Stainless steel saw blade for plaster of Paris.	-	WP-65	WP-65-2	WP-50	WP-50-2
Hardened steel saw blade for synthetic casts.	BS-77-2	BS-65	-	BS-50	BS-50-2
Hardened steel saw blade with Titan processing for synthetic casts.	-	BS-65-GT	-	BS-50-GT	-
Hardened steel saw blade with non-stick PTFE processing for synthetic casts.	-	BS-65-TF	BS-65-TF	-	-

VACUUM CLEANER FOR OSCIMED 2000

The vacuum cleaner is the same model like for the OSCIMED ERGO II, see page 5.

OSCIMED ERGO II

Cast cutter saw



ERGONOMIC, COMPACT, LIGHTWEIGHT AND EASY TO USE

OSCIMED ERGO II is a compact, lightweight plaster saw ergonomically shaped to enable precision work without any effort, even with only one hand. OSCIMED ERGO II is robustly built and equipped with a cooling fan for intensive use without overheating. The device is powered by a safe low voltage of 24 volts and complies with the latest medical standards.

Rated power	180 Watts
Length of the saw body	230 mm
Effective diameter	58 mm
Weight without power supply	860 gr.
Speed of motor rotation	14'500 rev/min
Noise level	67 dBA
Length of power supply cable	1,10 m
Length of the cable 24 VDC	3,00 m

OSCIMED ERGO Saw:

- with European plug: ref. OSC 340
- with British plug: ref. OSC 340-UK
- 120V/60 Hz Saw with U.S. plug: ref. OSC 340-US




Cooling fan for use without overheating.



A high-performance, easy-to-adapt dust vacuuming system.



Usable with one hand thanks to its very short (230 mm) and lightweight (860 g) shape.



Blade change without tools.

SAW BLADES AND ACCESSORIES



Blade lock nut
ref. OSC 202



Vacuum nozzle
ref. OSC 240

Saw Blades	Ø 77 mm	Ø 65 mm	Ø 65 mm	Ø 50 mm	Ø 50 mm
Stainless steel saw blade for plaster of Paris.	-	WP-65	WP-65-2	WP-50	WP-50-2
Hardened steel saw blade for synthetic casts.	BS-77-2	BS-65	-	BS-50	BS-50-2
Hardened steel saw blade with Titan processing for synthetic casts.	-	BS-65-GT	-	BS-50-GT	-
Hardened steel saw blade with non-stick PTFE processing for synthetic casts.	-	BS-65-TF	BS-65-TF	-	-

OSCIMED DUST EXTRACTION SYSTEM

4 levels of filtration, adjustable power and automatic start

Vacuum cleaner with adjustable suction power, mounted on four wheels. Its robust construction and power offer support for the saw.

The HEPA-compliant filters provide a standard filtration of 99.97% of particles of 0.3 microns. A power outlet allows the connection of the saw. The vacuum cleaner starts automatically upon activation of the saw and stops a few seconds after the saw stops.

Power: Adjustable from 300 to 1000 Watts

Noise level: 60-70 dB

Suction: 1800 mm water column, 40 litres/sec.

Length of suction hose: 2.5 meters.


Capacity of stainless steel tank: 20 litres.

Weight: 12 kg.

Filtration system with 4 levels

Standard HEPA:
99.97% of particles
of 0.3 microns

- ① HEPA filter ref. OSC 208-3
- ② HEPA filter ref. OSC 208-2
- ③ Polyester filter ref. OSC 208-1
- ④ Dust bag ref. OSC 207



Vacuum cleaner for OSCIMED ERGO II and OSCIMED 2000

OSCIMED vacuum cleaner with European plug: ref. OSC 206-V90

OSCIMED vacuum cleaner with British plug: ref. OSC 206-V90-UK

OSCIMED 120V/60 Hz vacuum cleaner with U.S. plug: ref. OSC 206-V90-US

ACCESSORIES



Vacuum cleaner hose
ref. OSC 272



Polyester filter
ref. OSC 208-1



Dust bag
ref. OSC 207



HEPA filter
ref. OSC 208-2



HEPA filter
ref. OSC 208-3



ERGONOMIC, LIGHTWEIGHT, COMPACT AND EASY TO USE

The OSCIMED PSV saw was designed for easy use, specifically for use with a vacuum cleaner. The rear mounting of the vacuum cleaner hose guarantees a perfect vision of the work area and ease of handling even for casts that are difficult to access. Its ergonomic design includes a grip area with reduced circumference. Its light weight allows single-handed use and the sensor switch is very easy to access.

Rated power	180 Watts
Length of the saw body	255 mm
Effective diameter	50 mm
Weight without power supply	890 gr.
Speed of motor rotation	14,500 rev/min
Noise level	63 dBA
Length of power supply cable	1.10 m
Length of the cable 24 VDC	3.00 m


Oscimed PSV Saw:

- with European plug: ref. OSC 300
- with British plug: ref. OSC 300-UK
- 120V/60 Hz Saw with U.S. plug: ref. OSC 300-US

MODERN TECHNOLOGY AND MAINTENANCE WITHOUT TOOLS

Its high-quality motor is based on vibration and noise dampers. The passage of the suctioned air ensures its cooling, so the saw can be used continuously without overheating. Accessories such as blades, vacuum nozzles and pipes can be changed rapidly, without tools. The device operates on a low voltage of 24 volts and meets the latest medical standards.

Change of blades without tools.



Simple sensor switch access.



Easy to use with one hand even with the dust extraction.

SAW BLADES AND ACCESSORIES



Blade lock nut
ref. OSC 202



Vacuum nozzle
ref. OSC 320

Saw Blades	Ø 77 mm	Ø 65 mm	Ø 65 mm	Ø 50 mm	Ø 50 mm
Stainless steel saw blade for plaster of Paris.	-	WP-65	WP-65-2	WP-50	WP-50-2
Hardened steel saw blade for synthetic casts.	BS-77-2	BS-65	-	BS-50	BS-50-2
Hardened steel saw blade with Titan processing for synthetic casts.	-	BS-65-GT	-	BS-50-GT	-
Hardened steel saw blade with non-stick PTFE processing for synthetic casts.	-	BS-65-TF	BS-65-TF	-	-

OSCIMED DUST EXTRACTION SYSTEM

4 levels of filtration, adjustable power and automatic start

Vacuum cleaner with adjustable suction power, mounted on four wheels. Its robust construction and power offer support for the saw.

The HEPA-compliant filters provide a standard filtration of 99.97% of particles of 0.3 microns. A power outlet allows the connection of the saw. The vacuum cleaner starts automatically upon activation of the saw and stops a few seconds after the saw stops.

Power: Adjustable from 300 to 1000 Watts

Noise level: 60-70 dB

Suction: 1800 mm water column, 40 litres/sec.

Length of suction hose: 2.5 meters.

Capacity of stainless steel tank: 20 litres.


Weight: 12 kg.



Filtration system with 4 levels

Standard HEPA:
99.97% of particles
of 0.3 microns

- ① HEPA filter ref. OSC 208-3
- ② HEPA filter ref. OSC 208-2
- ③ Polyester filter ref. OSC 208-1
- ④ Dust bag ref. OSC 207



OSCIMED vacuum cleaner for PSV saw with European plug: ref. OSC 206-V00
OSCIMED vacuum cleaner for PSV saw with British plug: ref. OSC 206-V00-UK
OSCIMED 120V/60 Hz vacuum cleaner for PSV saw with U.S. plug: ref. OSC 206-V00-US

ACCESSORIES



Vacuum cleaner hose
ref. OSC 311



Polyester filter
ref. OSC 208-1



Dust bag
ref. OSC 207



HEPA filter
ref. OSC 208-2



HEPA filter
ref. OSC 208-3

PLASTER TABLE

Thanks to its double articulation and to its low and wide frame, this plaster table is very sturdy.

Utilised together with the optional stand, this plaster table finds its place in every plaster room.

The dorsal part and the legs can be individually adjusted in every position thanks to the gas actuators.

The leg parts are adjustable separately till a down angle of 90°.

The table is adjustable in height thanks to an electrical actuator with foot or hand control.

The table is always equipped with castors and with a central locking which allows the castors to brake simultaneously.



ref. TP 01



Dimensions: 185 x 60 cm (without stand)
205 x 60 cm (with stand)

Adjustable height 65-115 cm

Capacity 175 kg

The frame is epoxy painted

Very little maintenance on the articulations and the mobile parts

European manufacturing

2 years warranty excepted on the upholstery

OPTIONS FOR PLASTER TABLE



Stand
ref. TP 06



Bar switch system for plaster table
ref. TP 04



Knee supports, included
adjustable side tubes
ref. TP 03



Additional charge for
the special width of 70 cm
instead of 60 cm
ref. TP 09

Extra footswitch
ref. TP 07



Extra handswitch
ref. TP 08

PLASTER TROLLEY



Position	Reference	Description
1	P-CHPL3P	Trolley fitted with: – 2 drawers – soaking compartment AVAILABLE COLOURS: sky blue or apple green
2	P-CHPL13P	Stand with Jersey distributor
3	P-PDCH3-167R	Movable workstation

PLASTER INSTRUMENTS



Cast spreader USA model, 23 cm
ref. MP-2010



Henning cast spreader, 28 cm
ref. MP-2034



Synthetic cast cutter,
serrated, 18 cm
ref. MP-1992-C



Stille cast shear, 26 cm
ref. MP-1997



Bruns bandage scissor,
plain bottom blade, 24 cm
ref. MP-1984



Lister bandage scissor
14 cm – ref. MP-1979
18 cm – ref. MP-1980



Bruns bandage scissor,
serrated bottom blade, 24 cm
ref. MP-1985



Synthetic cast cutter,
guarded end, serrated, 19 cm
ref. MP-1991-C

KNEE REST FOR CASTING

Knee rest with adjustable height
for an easy legs casting.




ref. OSC 355

STERILIZATION BASKETS

- High quality sterilization baskets in electropolished stainless steel AISI 304, size of the mesh 3x3x1 mm.
- double vertical frame for a high reliability of the wires welding.
- 2 handles.
- Stackable.
- No liquid stagnation!
- Lid with eccentric safety lock mechanism.






Basket complete with insert and fastening silicone bar.



Basket complete with lid and fastening silicone bar.

REFERENCES

	Model A	Model B	Model C
 Basket	Size 240x250x60 mm ref. STBK13-A	Size 480x250x60 mm ref. STBK13-B	Size 520x250x60 mm ref. STBK13-C
 Lid for basket with safety lock	Size 240x250 mm ref. STLD13-A	Size 480x250 mm ref. STLD13-B	Size 520x250 mm ref. STLD13-C
 Insert for basket with 2 handles	Size 205x230 mm ref. STIS13-A	Size 440x230 mm ref. STIS13-B	Size 480x230 mm ref. STIS13-C

FASTENING SILICONE BARS



IN



PUSH



CLIP



SLIDE

Model Clubs, length 210 mm for components of diameter 3-10 mm
ref. STSP-C



Clubs: Ø 3-10 mm - Blue

Model Snake, length 210 mm for components of diameter 0-10 mm
ref. STSP-S



Snake: Ø 0-3 mm - Blue



THE COMPANY OSCIMED SA

Oscimed was founded in year 2004. It is based in the swiss Jura, famous area well-known for its watch industries and high technologies. We focus on offering fine medical products for better daily living, at the best prices, with a high quality service, for the satisfaction of our customers.



OSCIMED SA – Allée du Quartz 3 – CH-2300 La Chaux-de-Fonds – Tel. +41 (0)32 926 63 36 – Fax +41 (0)32 926 63 50 – info@oscimedsa.com

La Chaux-de-Fonds

SWITZERLAND



Distributed by: