



Trioplus Advanced

A 3-section couch in the purest form - stable and functional. Extremely versatile: sitting, lying, lumbar flexion, Trendelenburg, manual traction, etc. Guarantees years of problem-free support for your manipulations.



- 3-section support surface, suitable for most basic positions.
- Adjustable middle section for the best possible position of the spinal column in flexion.
- Head support can be adjusted positively and negatively, and can also be adjusted almost fully downwards.
- Leg section also suitable for use as a back support.



Standard features

- Adjustable middle section.
- Central wheel-lifting mechanism with 4 casters.
- All sections gas spring supported.
- Gymna Stability Profile (GSP).
- Non slip feet.
- Cushions with padding of 40 mm.
- Breathing hole head section, incl. cushion.
- Cushions in a choice of 8 colours
- Prepared for hip traction support.
- Prepared for shoulder traction supports.



Specifications

Length x width	: ± 202 x 67cm
Height adjustment	: ± 45 - 95 cm
Weight	: ± 101 kg
Maximum Lifting capacity	: 150 kg
Electric motor strength	: 6.000 N
Hydraulic pump strength	: 10.000 N
Power supply	: 110 - 240 VAC 50/60 Hz or 100 VAC 50/60 Hz
Maximum consumption	: 1,5 A
Insulation class	: II
Moister resistance	: IP 54
Classification	: I
Conformity	: CE according to the Medical Devices Directive of the European Commission 93/42/EEG
Guarantee	: 3 Years (upholstery 1 year)

Ordering data - Trioplus Advanced



125 041
Electric Gymna
peripheral switch



125 042
Electric
footswitch



125 043
Hydraulic foot
pump

Accessoires

- 122 405 Paper roll holder, attaches to undercarriage
- 105 823 Paper rolls (per 12 pieces)
- 121 409 Extra footswitch for height adjustment
- 125 652 Hand switch for height adjustment
- 107 819 Traction support for hips or shoulders
- 120 040 Foot brace (lean back frame)
- 121 401 Gynaecological leg support
- 333 245 Traction device AutoTrac 460
- 121 434 Attachment frame for the AutoTrac 450/460
- 332 915 Foot end protector for cushion
Terry-cloth cover, stretch material (9 colours)

Optional

- 125 602 Electrical adjustable middle section
- 125 610 Motorized head section
- 121 404 Set brackets for fastening straps
- 121 402 Cushions 60 cm wide
- 121 403 Cushions 80 cm wide

