

# Reynard


## Disposable Underlays



### Disposable Underlays

Reynard Health Supplies 5 ply Underlays incorporate a non-woven facing, embossed polyethylene backing, a tissue absorbent layer and are sealed on all sides.

The underlays absorb urine and other body fluids while keeping the top layer dry. This keeps the users skin free from urine, body fluids and excremental material. The bottom sheet is completely impermeable and prevents side leakage.



#### Convenience

Available in 3 sizes

Large, Regular, Mini

#### Effective

Designed to resist shifting and bunching within the patients bed.

#### Cost-effective

Reynard's Disposable Underlays are cost—competitive.

#### Ideal for use in:

- Hospitals—Wards
- Hospital—Theatres
- Aged Care
- Private Hospitals
- Pharmacy
- Personal use
- Veterinary Care

**Directions:** Use the required number of Underlays  
Dispose of used underlays properly; do not flush.



Product specifications and availability may vary according to country.  
Always read the label and product information before use.

**Supporting documentation:** Full technical specifications are available upon request.  
**Bespoke solutions:** Should your market require more or fewer wipes, lightly scented alternatives or even biodegradable wipes, for example. Let us know.  
We are more than happy to develop bespoke solutions.

# Disposable Underlays




Reynard Disposable Underlays offer a highly absorbent material that holds moisture away from skin that is often sensitive. The absorbent layer contains different plies of tissue paper to assure high absorbency.



## RHS 503 Mini Disposable Underlays

### Technical specifications

- 30 x 20cm
- 50 Pieces per bag  
4 bags per carton



## RHS 505 Regular Disposable Underlays

### Technical specifications

- 43 x 60cm
- 50 Pieces per bag  
6bags per carton



## RHS 508 large Disposable Underlays

### Technical specifications

- 60 x 90cm
- 25 Pieces per bag  
4 bags per carton

- **Top Sheet:** soft and smooth, non-woven fabric—keeps skin dry and comfortable.
- **Absorbent layer:** Different piles of tissue paper assure high absorbency.
- **Bottom sheet:** Water proof PE film is absolutely impermeable and prevents side leakage.
- Guaranteed hygiene and cleanliness.

The recommendations regarding our preparations are based on scientific tests and are given in good faith. More detailed recommendations, e.g. Regarding material compatibility are only possible in particular cases. Our recommendations are without obligation and do not constitute a warranty. They do not preclude the company's own testing for the intended purposes and processes. In this respect we cannot accept any liability.

This complies with our general conditions of sale and supply.

### United Kingdom

Boston House  
Grove Technology Park  
Downsview Road, Wantage  
Oxfordshire OX12 9FF  
☎ +44 845 263 7335

### Australia

Suite 403, 44 Hampden Road  
Artarmon  
Sydney, NSW 2064  
enquiries @reynardhealth.com  
☎ 1300 667 199

### New Zealand

17 Napier Road  
Havelock North  
Hawke's Bay  
4130 New Zealand  
☎ +64 (0)6 650 0709

### International

17 Napier Road  
Havelock North  
Hawke's Bay  
4130 New Zealand  
☎ +64 (0)6 650 0709

[www.reynardhealth.com](http://www.reynardhealth.com)