


bed, more than 6 times what is required by the EN 60601-2-52 (2010) standard. You can be confident this bed will perform reliably throughout its useful life. With proper maintenance, it will remain safe for your patients and caregivers.

- Comparison performed on two different manufacturers' beds measuring the distance from the outer edge of the siderail to the frame in the upper position and the intermediate position
- 2. DJS Research, HR900 customer satisfaction survey, November 2011.
- To lower the rail, simply let it drop. You don't need to lower it manually. Multiplied by the number of times a caregiver can lower side rails, that can add up to over 200 kg of effort avoided each day.
- Markey, DW, Brown RJ, An interdisciplinary approach to addressing patient activity and mobility in the medical-surgical patient., Journal of Nursing Care Quality, 2002 Jul;16(4):1-12.
- Flanders SA, Harrington L, Fowler RJ. Falls and Patient Mobility in Critical Care: Keeping patients and staff safe. AACN Adv Crit Care. 2009; 20(3): 267-276.
- Cumulative load is the kg load multiplied by the number of cycles defined in each test protocol.

Tested to Hill-Rom, standards





 $Hill-Rom\ and\ Design\ and\ Hill-Rom\ word\ mark\ are\ trademarks\ or\ registered\ trademarks\ of\ Hill-Rom\ Services\ Inc.$ $Easy Chair\ is\ a\ trademark\ or\ registered\ trademark\ of\ Hill-Rom\ Services,\ Inc.$



Enhancing Outcomes for Patients and Their Caregivers....

France	+33 (0)2 97 50 92 12
United Kingdom	+44 (0)1530 411000
Deutschland	+49 (0)211 16450-0
Nederland	+31 (0)347 32 35 32
Italia	+ 39 (0)2 950541
Suisse/Schweiz	+41 (0)21 706 21 30
(deutschsprachig)	+41 (0)21 706 21 38


Österreich	+43 (0)2243 28550
Ireland	+353 (0)1 413 6005
Iberia	+34 (0)93 6856000
Nordic Countr	ies +46 (0)920 47 47 11
Export	Contact your local Hill-Rom
distributor or contact your Area Manager via	
website or call	+1 812 934 8173

www.hill-rom.com

Hill-Rom_® reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom_® makes is the express written warranty extended on the sale or rental of its products.

© 2012 Hill-Rom_® SARL, ALL RIGHTS RESERVED

Intuitive and efficient care: a better outcome.











on the side of the bed. The one-touch Egress position and the optimal side rail design allow for proper patient positioning the first time, every time.

less exertion for caregivers and greater confidence for both.



Bed low


With a single touch, you can use the Flat position function to move the bed from any articulated position to its lowest and flattest position – less than 39cm from the floor.



At a glance, you can be confident the bed is in the safest position for sleep – even from outside the patient's room.

The nightlight shines green when the bed is in the lowest position. It also illuminates the floor space where the patient exits the bed to help reduce the risk of falls.

If the bed is elevated from the lowest position, the nightlight will shine orange. A warning light on the control panel will also alert the caregiver.



engaged.



Safe sides

The side rails aid patients who are able to egress the bed independently. The Hill-Rom_o 900 bed is certified compliant to the EN 60601-2-52 hospital bed standard, meaning the risk of patient entrapment is minimised, even with the bed extension out.

