

EZ100



EZ Spray, catalog number EZ100, was designed for the rapid anesthetization of patients. Its single use, patented closed system prevents the possibility of cross-contamination and provides the benefit of deep, penetrating atomization.

EZ100 Features and Benefits

Low Cost Disposable System with Directional Tip and Calibrated Bottle


- Disposable, single use atomizer offers peace of mind by avoiding any possibility of cross-contamination
- Less expensive to use than processing and sterilizing a reusable system
- Directional tip provides health care workers with the ability to effortlessly direct and deliver medication to the desired location
- Calibrated, clear bottle allows for metered dose of medication. Bottles come in sizes ranging from 4 CCs to 100 CCs
- Fast and easy set up allows for immediate use

Patented Closed System

- Patented closed system adds another layer of protection to health care workers by preventing the aspiration and back flow of body fluid from patient and subsequent spread of possible pathogens
- Effortless delivery mechanism provides a controllable variation in the application rate of medication

Benefits of a Power Sprayer

- Incorporates a power sprayer system, which many doctors are already familiar with and like using
- The atomized medication deeply penetrates the nasal cavity, mouth, or throat
- It power sprayer features blows the atomized medication down below the cords
- Quick and thorough anesthetization
- Immediate anesthetization of patient allows for instant scoping
- Non-stick trigger assembly overcomes problems common to spraying crystallizing, viscous medications.



Ideals Areas of Use

- Anesthesia
- Pulmonary Function/Respiratory Care
- GI Lab
- Endoscopy
- Radiation Oncology



AlcoveMedical.com

800.709.9888