

seca sensa 804

The ability to shine,
without dazzling



- 24 electrodes for analysis of body fat and body water.
- Automatic user recognition for fast measuring results.
- Four memories for practical multi-person use.
- More than 170 years of experience for maximum precision.

seca sensa 804:

Body fat and body water analysis with new brilliance.

Take a look at the elegantly chromed electrodes – they are arranged to guarantee uninterrupted skin contact during the measuring of body fat and body water.

Chromed electrodes ensure a beautiful look.

A total of 24 electrodes show everything this personal scale can do: in addition to weight, it can also measure nutrition status. And it can do this extremely reliably since the electrodes' layout is not random, but very cleverly devised in the shape of a foot. Equally clever are the elegant design and the rounded-off weighing platform, which neither clashes with taste nor uncomfortable for your feet.

Finding time for those who are short on it.


Anyone wishing to measure BMI and/or body fat and body water can now do so in no time: simply enter personal data such as age, gender and fitness level once into one of the four memories – and the scale recognizes its users based on weight time and time again. If not, it asks the user by flashing the most likely memory on the display. Then simply select and you're finished.

You can't go wrong.

When designers create a seca scale, they consider not only the beautiful form, but also many practical details. In the case of the seca sensa 804 you'll find the rearfacing buttons which prevent data from being entered and results corrupted by accidental foot contact. You see, with a scale from seca nothing is left to chance.



The cleverly positioned buttons can't be accidentally activated while the scale is being used.



Chromed electrodes guarantee uninterrupted skin contact.

seca sensa 804

Technical data

- Capacity: 330 lbs / 150 kg
- Division: 0.2 lbs / 100 g
- Dimensions (WxHxD):
12.4 x 1.7 x 11.8" / 315 x 44 x 300 mm
- Weight: 2.9 lbs / 1.3 kg
- Functions: auto-HOLD, BMI, BF/BW-measurement, lbs/kg switch-over, user memory, automatic user recognition, tip-on automatic switch-on, automatic switch-off
- Power supply: batteries
- Parts to be delivered:
4 AAA batteries,
4 extension feet for carpets