

Sizing Guide

We recommend you take your actual body measurements because they are more accurate than over your clothes. Our garments have added stretch so that they will conform to your body.

You may wish to have someone help with this.


Hips: Standing feet together, measure around the fullest part of your bottom and over your hip bones.

Incontinence: The stretch in the undergarments will accommodate light to medium incontinence pads in your regular size. With heavy incontinence pads, we recommend measuring yourself with the incontinence pads in place for the most accurate sizing.


If you require a size outside our range please contact us to discuss.

DELLOCH HIPFIT STARTER KITS



2 HIPFIT PADS



3 PAIRS OF UNDERGARMENTS



MEASUREMENTS TAKEN AROUND WIDEST PART OF HIPS

	cm	inches	size
	79-85	31-33	XXS
	86-92	34-36	XS
	93-99	37-39	S
	100-106	40-41	M
	107-113	42-44	L
	114-120	45-47	XL
	88-94	35-37	XS
	95-101	38-40	S
	102-108	41-42	M
	109-115	43-45	L
	116-122	46-48	XL


FOR THE LATEST PRODUCT INFORMATION & TO BUY ONLINE
VISIT WWW.DELLOCH.COM



New Product Closed Pocket Protection

Ideal for dementia sufferers in facility care.

Closed Pocket Protection is ideal for people who want ready-to-go protection, no matter the time of day.



Comfortable in Bed	Sewn in Pads
Soft Foam Pad	Ready to Go
High Level of Protection†	Machine Washable



For more information please contact us on:

0800 658 814 | uslmedical.co.nz

uslcustomerservices@uslmedical.co.nz

Disclaimer: The Delloch Hip Protection System is designed to help prevent a fall related hip fracture. It is important that all the procedures for applying the pad in the correct position as outlined in the Quick Guide are followed in order to prevent a hip fracture as a result of a fall. However while this system can help lessen the severity of the impact of the fall, it does not guarantee that this will prevent a hip fracture.

Trademarks 2015 Delloch Limited
NZ Pat Appln No. 595893, NZ Pat Appln No. 595898,
NZ Design Reg No. 415437, NZ Design Reg No. 415438
The Delloch logo, Delloch Hip Protection, and other terms used in conjunction with this product are the trademarks of Delloch Ltd. Delloch is the owner of all past, present and future intellectual property used and developed in relation to Delloch Hip Protection and reserves all rights in relation to its intellectual property.



Prevent Hip Fracture Today

Stay Active & Independent



Redefining Excellence in Hip Protection

More Than Protection


Did you know that 90% of hip fractures are attributed to falls? Hip protectors have been proven to reduce the risk of hip fracture by nearly 60%*.

Delloch's advanced HipFit range provides you with the comfort and protection you need to step out into the world and continue the activities you love.

Our slimline protector pads have been designed to conform to your body whilst remaining discreet under clothes for invisible protection.

HipFit Adjustable Protection

A simple pad positioning system puts the protection in the right place.



Why Choose Delloch?

- ✓ **Outstanding Comfort** - Our cotton-rich fabrics and slimline pads provide everyday comfort.
- ✓ **Adjustable Protection** - Our unique HipFit System positions the pad in just the right place.
- ✓ **Ultimate Protection** - Our pads are thoroughly tested to ensure you receive all the protection you need.
- ✓ **Internationally Approved** - Our products are externally tested & approved for your safety.**
- ✓ **Choices for Your Lifestyle** - Our range allows you to choose the right product based on your needs.

*Source: Bentzen, H., Bergland, A., & Forsén, L. (2008). "Risk of hip fractures in soft protected, hard protected and unprotected falls" Inj. Prev., 14:306 – 310

**Tested to the SDMA Draft Standard: Hip Protectors Used As Medical Devices – Performance Requirements and Test Methods


Easy Wear & Easy Care

We believe these undergarments are the best you will own so we've gone the extra mile to deliver a product that not only looks great but feels great too.

- ✓ **Super soft cotton-rich fabrics** - Year-round comfort
- ✓ **Flat internal seams** - To prevent unwanted irritation
- ✓ **Soft waistbands** - For a gentle touch without rubbing
- ✓ **Flexible system** - All our pads fit all our undergarments


Female Classic Hips

For the women, the delicate detailing and tailored cut delivers an elegant, feminine style.



Male Classic Hips

For the men, the classic cut with channel knit ensures style as well as comfort.



Three Pad Options

Our pads fit all of our undergarments, so if your situation changes, so can your protection.

HipFit Plus

Ideal for people who are a medium to high falls or fracture risk.

HipFit Plus is ideal for people who need more protection or a little more reassurance.



Reinforced Polymer Shield	Comfort Curved Form
Soft Inner Foam Layer	Lightweight
Very High Level of Protection†	Easy to Wipe Clean

HipFit Active

Ideal for people who are a low to medium falls or fracture risk.

HipFit Active is a perfect match for those people wanting to stay active to retain health and mobility.



Comfortable In Bed	Comfort Curved Form
Dual Density Foam Pad	Ultra Lightweight
High Level of Protection†	Easy to Wipe Clean

HipFit Flexi

Ideal for people who wish to wear hip protectors in bed.

HipFit Flexi has been specifically developed to be comfortable when lying on and is a great choice at night.



Comfortable In Bed	Machine Washable
Soft Foam Pad	Ultra Lightweight
High Level of Protection†	Easy to Wipe Clean

†Tested in accordance to the International Hip Protection Research Group: Recommendations for Bio-Mechanical Testing of Hip Protectors.