



Welch Allyn  
SureTemp® Plus 690 Electronic Thermometer



Taking temperatures is easier, safer and more accurate—so taking care of all your patients is more convenient.

# The SureTemp Plus 690—with the speed and accuracy you require and the versatility you need for use with all your patients, all the time.

## 1 Easy:

An easy-to-read LCD display with flashing icons that clearly communicate in any language, and an intuitive user interface and ergonomic shape

## 2 Safe:

Patient safety is improved thanks to an interchangeable probe well that helps reduce risk of cross-contamination

## 3 Fast and Accurate:

In addition to improved accuracy, the SureTemp Plus 690 is one of the fastest thermometers on the market, providing approximately 4-6 second oral readings, 10-second paediatric axillary and rectal readings and 15-second adult axillary readings

## 4 Versatile:

To top it all off, you can now use the SureTemp Plus 690 with all your patients, from infants to adults—whether you need to take oral, axillary or rectal temperatures

- Easy-to-use, intuitive operation
- Interchangeable oral/axillary and rectal probe wells help reduce the risk of cross-contamination



# The **SureTemp Plus 690**

Removable probe well for easy cleaning

User-selectable icon and/or words for identification of thermometer mode

C/F button for toggling between Fahrenheit and Celsius

Prominent recall button for display of last temperature taken

Ergonomic design with rubberised grips

Convenient storage for 25 probe covers

Large, easy-to-read LCD

Battery life indicator

Mode button for selecting oral, paediatric axillary or adult axillary



The **SureTemp Plus 690** with optional wall mount has storage space for 50 additional probe covers.

# The SureTemp Plus 690—setting the standard of care even higher

SureTemp Plus 690 Product Specifications	
Calibration Accuracy	±0.1 °C Monitor Mode
Temperature Predict Time	Oral: Approximately 4-6 seconds Adult Axillary (18 years and older): Approximately 12 - 15 seconds Paediatric Axillary (17 years and younger): Approximately 10 - 13 seconds Rectal: Approximately 10 - 13 seconds
Unit Dimensions	215 mm x 81 mm x 62 mm
Weight	357 g
Ambient Temperature Range	10 °C to 40.0 °C
Patient Temperature Range	26.7 °C to 43.3 °C
Power Source	(3) 1.5 Vdc AA batteries
Battery Life	Approximately 6,000 temperatures

## Ordering Information:

### THERMOMETERS

- 01690-401** SureTemp Plus 690 Electronic Thermometer, Oral Probe and Oral Probe Well with 1.2 m cord  
**01690-300** SureTemp Plus 690 Electronic Thermometer, One per room, Wall Mount, Oral Probe and Oral Probe Well with 2.7 m cord (must be locked to the wall)

### ACCESSORIES

- 21326-0000** Wall Holder  
**21330-0000** Holder for Spare Probe and Well  
**06138-0000** Calibration Key Assembly 690/692  
**02892-000** Probe and Well Kit, 1.2 m Rectal  
**02892-100** Probe and Well Kit, 2.7 m Rectal  
**02893-000** Probe and Well Kit, 1.2 m Oral  
**02893-100** Probe and Well Kit, 2.7 m Oral  
**02891-0000** Probe Well Blue  
**02891-1000** Probe Well Red
- 08273-000** Thermometer Stand with Holder for Spare Probe and Well  
**21333-0000** Bracket, M600 Stand Upgrade with Holder for Spare Probe and Well

### REPLACEMENT PROBE COVERS

- 05031-101** Case of 1,000 Probe Covers  
**05031-105** Case of 5,000 Probe Covers

To find out more, please contact your **Welch Allyn Sales Representative** or Customer Service on **1800 650 083**



Welch Allyn Australia Pty. Ltd.  
 Suite 4.01, 2-4 Lyonpark Road  
 Macquarie Park NSW 2113 AUSTRALIA