

Flocare® Feeding Sets

Flocare® Feeding Sets are final ENFit™ and designed to safely and securely connect feeds to feeding tubes and accessories

Features

- Feeding sets will connect to ENFit™ syringes at the medication port for increased patient safety.
- Flocare® pack feeding sets connect to Nutrison and Nutrini packs by a plus shaped spike end (ENPlus), which can only be connected to enteral feed packs and containers, increasing patient safety.
- Flocare® universal feeding sets connect to Nutrison and Nutrini bottles by a screw connector.
- Flocare® Mobile Giving Sets are available with a convenient, shorter length of tubing to support mobile feeding.
- All Flocare® feeding sets are compatible with the Flocare® Containers (40441, 40455, 70139)
- Latex and DEHP free materials
- Coloured lilac to make a clear distinction from parenteral feeding systems

Standard Giving Sets (30 Units per Carton)

For use with ENPlus tube feed packs:



86514
Infinity Pack
Set 3-way
port & drip
chamber



95347
Infinity Pack
Set Y-port



95349
Infinity Pack Set
no port



90084
Infinity Bottle Set
3-way port & drip
chamber



101587
Infinity Bottle Set
Y-port & drip
chamber



90086
Infinity Bottle Set
no port & drip
chamber

Mobile Giving Sets (30 Units per Carton)

For use with ENPlus tube feed packs:



101586
Infinity Pack
Mobile Set Y-
port



86484
Infinity Pack
Mobile Set no
port



86460
Gravity Pack Set
Y-port & drip
chamber



89755
Gravity Pack &
Bottle Set
Y-port & drip
chamber

Gravity Giving Sets (30 Units per Carton)

Not for use with tube feeding pumps

Recommendations for Use

- Flocare® feeding sets are intended for 'single use only' and should be changed every 24 hours

Ordering Information

- To order contact
Nutricia Customer
Care 1800 889 480

- ✓ Enhancing patient safety
- ✓ Increasing simplicity
- ✓ Reducing risk of misconnections



COMPLETE RANGE, COMPLETE SUPPORT, COMPLETE SAFETY

FLOCARE Giving Sets - Transition		FLOCARE ENFit Giving Sets - Final	
Current Code	New Code	New Product Description	Transition Status Australia
Standard Giving Sets for Flocare Infinity and Infinity ⁺ Pumps			
For use with ENPlus tube feed packs			
91403	86514	Flocare Infinity Pack Set - 3-way port & drip chamber	✓ Completed
95168	95347	Flocare Infinity Pack Set - Y port	✓ Completed
95174	95349	Flocare Infinity Pack Set - no port	✓ Completed
For use with enteral bottles			
91408	90084	Flocare Infinity Bottle Set - 3-way port & drip chamber	✓ Completed
95172	101587	Flocare Infinity Bottle Set - Y port & drip chamber	✓ Completed
95179	90086	Flocare Infinity Bottle Set - no port & drip chamber	✓ Completed
Mobile Giving Sets for Flocare Infinity and Infinity ⁺ Pumps			
95181	101586	Flocare Infinity Pack Mobile Set - Y port	✓ Completed
95182	86484	Flocare Infinity Pack Mobile Set - no port	✓ Completed
Gravity Sets – Not for use with enteral feeding pumps			
91402	86460	Flocare Gravity Pack Set - Y port & drip chamber	✓ Completed
91405	89755	Flocare Gravity Pack & Bottle Set - Y port & drip chamber	✓ Completed

NEW ENFit FLOCARE FEEDING ACCESSORIES &SYRINGES

Product	Code
ENFit Transition Adaptor	122073
Flocare 5 Step Adaptors 5 x 6 packs	89828
Flocare ENFit Bolus Adaptor	89740
Flocare ENFit Feeding Connector Ch18	91396
Infinity Drip Stand	115591
Flocare Back Pack (Adult)	115603
Flocare Back Pack (Child)	115592
ENFit Syringe 10ml (Single Use)	118498
ENFit Syringe 20ml (Single Use)	118497
ENFit Syringe 60ml (Single Use)	118496
ENFit Syringe 10ml (7 Day Use)	119408
ENFit Syringe 20ml (7 Day Use)	119409
ENFit Syringe 60ml (7 Day Use)	119410

NEW ENFit FLOCARE NASO-GASTRICTUBES

Product	Code
Flocare PUR ENFit NG Tube Ch5 50cm	94566
Flocare PUR ENFit NG Tube Ch8 110cm	94569
Flocare PUR ENFit NG Tube Ch10 110cm	94570
Flocare PUR ENFit NG Tube Ch12-110cm	94571

NEW ENFit FLOCARE Balloon GASTROSTOMYTUBES

Product	Code
Flocare ENFit G tube Ch14	94815
Flocare ENFit G tube Ch18	94817
Flocare ENFit G tube Ch20	94818